CZ - INFORMACE PRO UŽIVATELE

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EN - USER INFORMATION

FLOWER - powder free nitrile examination gloves blue, non-sterile

Available size: S, M, L, XL


1) Medical Device Regulation (MDR) (EU) 2017/745

- This product is classified under Class I Medical Device per Rule 1 and Rule 5 of Annex IX, meets the provisions of the Medical Device Regulation (MDR) (EU) 2017/745.
- This product complies with European Standards EN 455-1:2000, EN 455-2:2015, EN 455-3:2015, and EN 455-4:2009.

2) PPE EU Type-Examination

- This product is classed as Category III Personal Protective Equipment (PPE) according to PPE Regulation (EU) 2016/425 and has been shown to comply with this Regulation through the Harmonised European Standards EN 420:2003+A1:2009, EN ISO 374-1:2016/Amd.1:2018(E), EN ISO 374-1:2016 and EN ISO 374-5:2016.
- Notified Body responsible for certification is SATRA Technology Europe Limited (2777), Bracetown Business Park, Clonee, Dublin 15, D15 YN2P, Ireland.
- Notify Body responsible for internal production control plus supervised product checks at random intervals (Module C2) is SATRA Technology Europe Limited (2777), Bracetown Business Park, Clonee, Dublin 15, D15 YN2P, Ireland.
- The EU Declaration of Conformity is accessible at www.batist.com

3) Marking

- Micro Organism Hazards Pictogram: EN ISO 374-5:2016 Protect against Bacteria, Fungi and Virus. No penetration of bacteriophages through the specimen and the following pictogram is applied.
- Chemical Hazards Pictogram: EN 16523-1:2015; Additional information on chemical resistance obtainable from manufacturer.

EN ISO 374-1:2016 permeation levels are based on breakthrough times as follows:

| Permeation Performance Level | 1 | 2 | 3 | 4 | 5 | 6 |
|----------------------------------|-----|-----|-----|------|------|------|
| Measured breakthrough time (min) | >10 | >30 | >60 | >120 | >240 | >480 |

This product complies with Type B requirements and the following pictogram shall be used with reference to clause 6.2 of ISO 374-1.

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- The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.
- The gloves were found to meet with the REACH annex XVII requirements for Polycyclic Aromatic Hydrocarbons (PAHs).
- Components used in glove manufacturing may cause allergic reactions in some users. If allergic reactions occur, seek for medical advice immediately.

5) Gloves for Special Applications (EN 420:2003, Clause 5.1.3)

These gloves are designed to protect the hand surface from chemical splashes. Therefore, the length of the gloves is below EN requirements of total minimum glove length, and deems as 'Fit for Special Purpose'.

6) Product Instruction for Use

- Usage** - For Single Use only. If re-used, the risk of contamination and infection increases due to improper cleaning processes; and increased risk of holes and tear during re-use due to weakening of gloves by cleaning processes.
- Sizing** - Select the right size glove for your hand.
- Donning** - Hold glove by the bead with one hand. Align the glove thumb with your other hand thumb and slide your hand into the glove, one finger into each glove finger. Pull by the glove palm to a get a good fit. Don the other glove by the same procedure.
- Inspection** - Punctures or tears may occur after donning. Inspect each glove after donning, and immediately discontinue use if found damaged.
- Doffing** - Hold glove bead and pull toward the finger until the glove come off.
- Disposal** - Properly disposal of all used gloves. Follow your Institution's policies for disposal.
- Handling and Storage**
Store in a cool and dry place. Opened boxes should be kept away from fluorescent and sunlight. Gloves are packed in dispenser which is suitable for transport. Keep the gloves in the box when not in use.
- Shelf life**
The shelf life of product is 3 years from date of manufacture.

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