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[Use for] Mask is made from nonwoven textile meltblown and filtration layer. Nose wire is from iron coated with polypropylene. 2 earloops for fixation.

[Size] 16 x 10,8 cm

[Standard] EN 149:2001+A1:2009

[Indications For Use]

The protective mask is intended used for protection against solid and oily particulates, liquid and microorganisms such as coal dust, cement dust, acid fog, paint fog, oil smoke, oil mist, asphalt smoke, coke oven smoke etc.
Single-use.

[Fitting]

- Pull an extend the mask inner side to the face with two hands ensure the nose clip is upward.
- Hold the mask on the face and pull the earloop to the ears.
- Press the nose clip to fit the shape of the nose.
- Press the mask to fit the respirator onto your face.

[USE CAUTIONS]

- Please check whether the packaging is damaged, whether the mask body is damaged, deformed or has other obvious defects before use.
- It is unlikely that the requirements for leakage will be achieved if facial hair passes under the face seal.
- Please do not use the mask under low air quality condition, such as serious contaminants in air, oxygen deficiency.
- Do not use the mask in explosive atmosphere;
- After use, please discard the mask following local regulation.
- The mask is only for single use, do not use repeatedly.

Read instructions for use

<60% Store in dry place

Non reusable

-20C 30C Temperature range of storage conditions

Personal protective equipment category III

FFP2 – protection class
NR – single use

Longyan Branch of Kangyuan Health Technology Co., Ltd.
Building 9, No.1 Tengfei 6th Road, ChangTing Country,
Fujian, China.

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Luxus Lebenswelt GmbH,