CZ - INFORMACE PRO UŽIVATELE

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EN - USER INFORMATION

IMMUNITY – LATEX powder-free surgical gloves, sterile

Available size: 5, 5, 6, 6, 5, 7, 7, 5, 8, 8, 5, 9

1) MD - Medical Device Directive

- This product is classified under Class IIa Medical Device per Annex IX, meets the provisions of the Medical Device EC Council directive 93/42/EEC.
- Notified Body responsible for certification is TÜV SÜD Product Service GmbH (0123), Ridlerstrasse 65, 80339 München, Germany
- This product complies with European Standards EN 455-1:2000, EN 455-2:2015, EN 455-3:2015, and EN 455-4:2009.

2) PPE EU Type-Examination

- This product is classified as Category III Personal Protective Equipment (PPE) according to PPE Regulation (EU) 2016/425 and has been shown to comply with this Regulation through the Harmonised European Standards EN 420:2003+A1:2009, EN ISO 374-1:2016, EN ISO 374-2:2014, EN ISO 374-4:2013 and EN ISO 374-5:2016.
- Notified Body responsible for certification/PPE is SATRA Technology Europe Limited (2777), Bracetown Business Park, Clonee, Dublin 15, D15 YN2P, Ireland.
- Notify Body responsible for internal production control plus supervised product checks at random intervals (Module C2) is SATRA Technology Europe Limited (2777), Bracetown Business Park, Clonee, Dublin 15, D15 YN2P, Ireland.
- The EU Declaration of Conformity is accessible at www.batist.com

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
3) Marking

- Micro Organism Hazards Pictogram: EN ISO 374-5:2016 Protect against Bacteria, Fungi and Virus. No penetration of bacteriophages through the specimen and the following pictogram is applied.

- This product complies with Type B requirements and the following pictogram shall be used with reference to clause 6.2 of EN ISO 374-1:2016 .

4) Performance and Limitation of Use

- This product had been tested in accordance with EN ISO 374-5:2016 to provide protection against Bacteria, Fungi and Virus with result: Protection against bacteria and fungi - Pass Protection against viruses - Pass
- The gloves had been tested in accordance with EN 16523-1:2015 resistance to permeation by chemicals and achieved the following performance levels:

Chemicals	Performance Level
40% Sodium Hydroxide (K)	6
37% Formaldehyde (T)	6
30% Hydrogen peroxide (P)	6

EN ISO 374-1:2016 permeation levels are based on breakthrough times as follows:

Permeation Performance Level	1	2	3	4	5	6
Measured breakthrough time (min)	>10	>30	>60	>120	>240	>480

- This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.
- The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm² – where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical used in a mixture.
- It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.
- When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.
- Before usage, inspect the gloves for any defect or imperfections.
- This product had been tested in accordance with EN 374-4:2013 and achieved the following degradation results:

Chemicals	Mean Degradation / %
40% Sodium Hydroxide (K)	-11
37% Formaldehyde (T)	2,6
30% Hydrogen peroxide (P)	10,1

- EN 374-4:2013 Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemicals.
- This product provides protection against Bacteria, Fungi and Virus.
- The gloves had been tested in accordance with ISO 16604:2004 to meet the requirements of EN ISO 374-5:2016 for resistance to penetration by blood-borne pathogens-test method using Phi-X174 bacteriophage.
- The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.
- The gloves were found to meet with the REACH Annex XVII requirements for Polycyclic Aromatic Hydrocarbons (PAHs).
- The gloves contain natural rubber latex, which may cause allergic reactions or anaphylactic shock to sensitive individuals. If you have experience with reaction to latex, don't use the product and immediately consult your physician. Also in case of occurrence irritation or allergic reactions, contact medical emergency immediately.**

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5) Product Instruction for Use

- Usage** – For Single Use only. If re-used, the risk of contamination and infection increases due to improper cleaning processes; and increased risk of holes and tear during re-use due to weakening of gloves by cleaning processes.
 - Sizing** – Select the right size glove for your hand.
 - Donning** – Hold glove by the bead with one hand. Align the glove thumb with your other hand thumb and slide your hand into the glove, one finger into each glove finger. Pull by the glove palm to a get a good fit. Don the other glove by the same procedure.
 - Inspection** – Punctures or tears may occur after donning. Inspect each glove after donning, and immediately discontinue use if found damaged.
 - Doffing** – Hold glove bead and pull toward the finger until the glove come off.
 - Disposal** – Properly disposal of all used gloves. Follow your Institution's policies for disposal.
- #### 6) Handling and Storage
- Store in a cool and dry place. Keep away from fluorescent and direct sunlight. Keep the gloves in the original packaging when not in use. Open packaging just before use.
- #### 7) Shelf life
- The shelf life of product is 5 years from date of manufacture.

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