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**EN - USER INFORMATION****FLOWER – LATEX powder-free examination gloves, non-sterile**

Available size: S, M, L, XL

**1) Medical Device Regulation (MDR) (EU) 2017/745**


- This product is classified under Class I Medical Device per Rule 1 and Rule 5 of Annex IX, meets the provisions of the Medical Device Regulation (MDR) (EU) 2017/745.
- This product complies with European Standards EN 455-1:2000, EN 455-2:2015, EN 455-3:2015, and EN 455-4:2009.

**2) PPE EU Type-Examination**

- This product is classed as Category III Personal Protective Equipment (PPE) according to PPE Regulation (EU) 2016/425 and has been shown to comply with this Regulation through the Harmonised European Standards EN 420:2003+A1:2009, EN ISO 374-1,2,4,5
- Notified Body responsible for certification is SATRA Technology Europe Limited (2777), Bracetown Business Park, Clonee, Dublin 15, D15 YN2P, Ireland.
- Notified Body responsible for quality assurance of the production process (Module D) is SGS FIMKO OY (0598 P.O. Box 30 (Särkinientie 3), 00211 Helsinki, Finland
- The EU Declaration of Conformity is accessible at [www.batist.com](http://www.batist.com)


**3) Marking**

- Micro Organism Hazards Pictogram: EN ISO 374-5:2016 Protect against Bacteria, Fungi and Virus. No penetration of bacteriophages through the specimen and the following pictogram is applied.

VIRUS

- This product complies with Type C requirements and the following pictogram shall be used with reference to clause 6.3 of ISO 374-1.

ISO 374-1 / Type C

**4) Performance and Limitation of Use**

- This product has been tested in accordance with EN ISO 374-5:2016. Protection against bacteria and fungi - Pass Protection against viruses - Pass
- Gloves had been tested in accordance with EN 16523-1:2015 resistance to permeation by chemicals and achieved the following performance levels:

| Chemicals                | Performance Level |
|--------------------------|-------------------|
| 40% Sodium Hydroxide (K) | 3                 |
| 37% Formaldehyde (T)     | 5                 |
| 96% Sulphuric acid (L)   | 1                 |

EN ISO 374-1:2016 permeation levels are based on breakthrough times as follows:

| Permeation Performance Level     | 1   | 2   | 3   | 4    | 5    | 6    |
|----------------------------------|-----|-----|-----|------|------|------|
| Measured breakthrough time (min) | >10 | >30 | >60 | >120 | >240 | >480 |

- This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.
- The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm – where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical used is a mixture.
- It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.
- When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.
- Vf usage, inspect the gloves for any defect or imperfections.
- This product had been tested in accordance with EN 374-4:2013 and achieved the following degradation results:

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- EN 374-4:2013 Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemicals.
- This product provides protection against Bacteria, Fungi and Virus.
- The gloves had been tested in accordance with ISO 16604:2004 to meet the requirements of EN ISO 374-5:2016 for resistance to penetration by blood-borne pathogens-test method using Phi-X174 bacteriophage.
  - The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.
- The gloves contain natural rubber latex, which may cause allergic reactions or anaphylactic shock to sensitive individuals. If you have experience with reaction to natural latex, don't use the product and immediately consult your physician. Also in case of occurrence irritation or allergic reactions, contact medical emergency immediately.**

**5) Gloves for Special Applications (EN 420:2003+A1:2009)**

These gloves are designed to protect the hand surface from chemical splashes. Therefore, the length of the gloves is below EN requirements of total minimum glove length, and deems as 'Fit for Special Purpose'.

**6) Product instruction for use****a) Rules for glove use in the medical field:**

- Wear the correct gloves when needed.
- Wear gloves no longer than 2 hours.
- Wash hands once gloves have been removed.
- Disposable gloves must be discarded once removed. Do not save for future use.
- Remove gloves before touching personal items, such as phones, computers, pens and one's skin.
- Do not wear gloves out of working office. If gloves are needed to transport anything, wear one glove to handle to transported item.
- If for any reason a glove fails, and chemicals come into contact with skin, consider it an exposure and seek medical attention.

**b) Find the right size** - Find the appropriate size of the gloves. User can determine if their gloves too small by extending their fingers until straight.**c) Donning:**

- Remove all hand and wrist jewelry, and wash the hands before donning.
- Place the gloves on the prepared work surface.
- The user puts a glove on his/her dominant hand by grabbing it with the other hand, remembering to only touch the inside of the gloves, and slipping it over the dominant hand until it reaches finger level.
- The wearer uses the gloved dominant hand to slip the other glove onto the non-dominant hand.
- Once both glove are on, the users can touch the outside of the gloves to ensure a proper fit.

**d) Doffing:**

- Using the dominant hand, users start by grabbing the outside of the glove on the non-dominant hand on the palm side near the cuff.
- Pull the glove off the non-dominant hand and place it in the gloved hand, balling it up.
- Slip two fingers under the cuff of the other hand glove and carefully peel it off the hand without touching the wrist, turning the remaining glove inside out as it is removed and in turn encasing the first glove.

**e) Disposal** - Properly disposal of all used gloves. Follow your Institution's policies for disposal.**7) Handling and Storage**

Store in a cool, dry and well-ventilated area away from extremes of temperature. Opened boxes should be kept away from fluorescent and sunlight. Gloves are packed in dispenser which is suitable for transport. Keep the gloves in the box when not in use.

**8) Shelf life**

The shelf life of product is 5 years from date of manufacture.

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