Paper size : 420mm x 297mm (A3)
Fold size : 210mm x 297mm (A4)

Fold line

CZ – INFORMACE PRO UŽIVATELE

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EN – USER INFORMATION

FLOWER – LATEX powder-free examination gloves, non-sterile
Available size: XS, S, M, L, XL

1) Medical Device Regulation (MDR) (EU) 2017/745

- a) This product is classified under Class I Medical Device per Rule 5 of Annex VIII, meets the provisions of the Medical Device Regulation (MDR) (EU) 2017/745.
- b) This product complies with European Standards EN 455-1:2000, EN 455-2:2015, EN 455-3:2015, and EN 455-4:2009.

2) PPE EU Type-Examination

- a) This product is classed as Category III Personal Protective Equipment (PPE) according to PPE Regulation (EU) 2016/425 and has been shown to comply with this Regulation through the Harmonised European Standards EN ISO 21420:2020, EN ISO 374-1:2,4,5
- b) EU Type Examination by SATRA Technology Europe Limited (NB 2777) Braetown Business Park, Clonee, D15 YN29 Ireland.
- c) Conformity assessment based on Regulation (EU) 2016/425 Annex VIII (Module D) by SGS FMKO OY (NB 0598) Takomite 6, FK00380, Heikinki, Finland
- d) The EU Declaration of Conformity is accessible at www.batist.com

3) Marking

- a) Micro Organism Hazards Pictogram: EN ISO 374-5:2016 Protect against Bacteria, Fungi and Virus. No penetration of bacteriophages through the specimen and the following pictogram is applied.

EN ISO 374-5:2016

EN ISO 374-5:2016

- b) This product complies with Type C requirements and the following pictogram shall be used with reference to clause 6.3 of EN ISO 374-1.

EN ISO 374-1:2016-A1:2018 / Type C

4) Performance and Limitation of Use

- a) This product has been tested in accordance with EN ISO 374-5:2016. Protection against bacteria and fungi - Pass
- b) Protection against viruses - Pass
- c) Gloves had been tested in accordance with EN 16523-1:2015 resistance to permeation by chemicals and achieved the following performance levels:

Chemicals	Performance Level
40% Sodium Hydroxide (K)	3
37% Formaldehyde (T)	5
96% Sulphuric acid (L)	1

EN ISO 374-1:2016-A1:2018 permeation levels are based on breakthrough times as follows:

Permeation Performance Level	1	2	3	4	5	6
Measured breakthrough time (min)	>10	>30	>60	>120	>240	>480

- i) This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.
- ii) The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the cuff is equal to or over 400 mm² - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical used is a mixture.
- iii) It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.
- iv) When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.
- v) Before usage, inspect the gloves for any defect or imperfections.
- c) This product had been tested in accordance with EN ISO 374-4:2019 and achieved the following degradation results:

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