

INFECTION PREVENTION

Surface disinfection – Wet wipes



DISINFECTANT WIPES



Effective
Against
Coronavirus

Disinfectant Wipes offer great flexibility. Powerful combination chemistry delivers broad spectrum efficacy that's ideal for surfaces and non-invasive medical devices.

- Ideal for all surfaces made of materials that are not resistant to alcohol.
- Used for simultaneous cleaning and disinfection of surfaces.
- Ready to use wipes for disinfection surface of noninvasive medical devices.
- Alcohol free wipes.
- For disinfecting all touch surfaces.
- Tested in accordance with EN 16615 – “Four Field Test”
- Microbiological activity from 1 minute against bacteria, yeast and fungi, mycobacteria and viruses (including coronavirus)

Active Ingredients:

N,N-Didecyl-N,N-Dimethylammonium Chloride
N-(3Aminopropyl)-N-dodecylpropane-1,3diamine

Directions for use:

Use one wipe to disinfect the Surface. Allow to dry naturally. On heavily soiled surfaces, more than one wipe may be required. Do not re-use wipes. Observe the duration of action to achieve the microbiological effect.

Packing:



REF: 1320100497

Product code: S524110MPCE

Medipal Disinfectant Wipes 50 Flow-Wrap
(50 wipes in Flow-Wrap). 24 pcs in case.
Size of single wipe: 200 x 175 mm



REF: 1320100498

Product code: S522110MPCE

Medipal Disinfectant Wipes 200 Flow-Wrap
(200 wipes in Flow-Wrap). 6 pcs in case.
Size of single wipe: 250 x 200 mm

Shelf life 3 years. After opening 1 month.

INFECTION PREVENTION

Surface disinfection – Wet wipes



ALCOHOL WIPES

These highly efficacious general purpose wipes Isopropylalcohol (IPA). With an ultra-low linting material which delivers a controlled dosage to ensure optimum disinfection performance with every wipe.

- Ready to use wipes for disinfection surface of noninvasive medical devices.
- Ideal for all surfaces made of materials that are resistant to alcohol.
- Contains 70% isopropyl alcohol
- They do not leave streaks on surfaces
- For disinfecting all touch surfaces and small noninvasive medical equipment
- Tested in accordance with EN 16615 – “Four Field Test”
- Microbiological activity in 60 seconds against bacteria, yeast, mycobacteria and viruses (including Noro viruses)

Active Ingredients:

Isopropanol 70% v/v

Directions for use:

Use one wipe to disinfect the Surface. Allow to dry naturally. On heavily soiled surfaces, more than one wipe may be required. Do not use on surfaces sensitive on alcohol . Do not re-use wipes. Do not use on patient and on invasive equipment. After use allow the surface to dry naturally before your next procedure. Observe the duration of action to achieve the microbiological effect.

Packing:



REF: 1320100495

Product code: W600110MPCE

Medipal Disinfectant Wipes 200 Tub
(200 wipes in Tub). 10 pcs in case.

Size of single wipe: 195 x 190 mm

Shelf life 3 years. After opening 6 months.

INFECTION PREVENTION

Surface disinfection – Wet wipes



CHLORINE WIPES



Our ready-to-use Chlorine Wipes are less toxic and less corrosive than other chlorine-releasing products whilst delivering a rapid biocidal action against a broad spectrum of micro-organisms and spores across non-porous hard surfaces, equipment and non-invasive medical devices.

- Ready-to-use wipes for disinfection of non-porous hard surfaces, equipment and non-invasive medical devices.
- Rapid biocidal action against a broad spectrum of micro-organisms and spores
- Tested in accordance with EN 16615 – “Four Field Test”
- Microbiological activity in 60 seconds against bacteria, mycobacteria and viruses (including coronaviruses).
- Full range of microbiological activity including spores in 2 minutes

Active Ingredients:

Sodium Hypochlorite

Directions for use:

Wear protective gloves. Use one wipe to disinfect the surface. Allow to dry naturally. On heavily soiled surfaces, more than one wipe may be required. Do not re-use wipes. After use, allow the surface to dry naturally before your next procedure. Observe the duration of action to achieve the microbiological effect.

Packing:



REF: 1320100496

Product code: S570110MPCE

Medipal Chlorine Disinfectant Wipes 50 Flow-Wrap
(50 wipes in Flow-Wrap). 12 pcs in case.
Size of single wipe: 200 x 200 mm

Shelf life 1 year. After opening 1 month.

INFECTION PREVENTION

Surface disinfection – Wet wipes



Microorganism	NONALCOHOL WIPES		ALCOHOL WIPES
	Medipal Disinfectant Wipes	Medipal Chlorine Disinfectant Wipes	Medipal 70% IPA Wipes
Bactericidal Efficacy	Contact time		
EN 13727	30 s	1 min	1 min
EN 13697	1 min	2 min	30 s
EN 16615	1 min	1 min	1 min
Yesticidal Efficacy			
EN 13624	1 min	2 min	1 min
EN 13697	1 min	2 min	
EN 16615	1 min	1 min	1 min
EN 14562	2 min		
Fungicidal Efficacy			
EN 13624		2 min	
Mycobactericidal Efficacy			
EN 14348	1 min	1 min	1 min
EN 14563	1 min		1 min
Virucidal Efficacy			
EN 14476			
Enveloped viruses	1 min	1 min	5 min
Limited virucidal	1 min	1 min	1 min*
Full virucidal	2 min	1 min	
Sporocidal Efficacy			
EN 13704			
Clostridium difficile		2 min	
Clostridium sporogenes		2 min	
Bacillus subtilis		2 min	
Bacillus cereus		2 min	

*Murine Norovirus

How to use:

- STEP 1** Please follow your agreed risk assessment policy guidelines regarding the use of PPE.
- STEP 2** Choose the dispenser type and dispense the wipe.
- STEP 3** Wipe the surface in an S Shape moving from clean to dirty. Use the wipe flat not scrunched. Do not go over the same area twice with the same wipe.
- STEP 4** Use a fresh wipe if your wipe becomes soiled or dry.
- STEP 5** Discard used wipes in the appropriate waste bin following your local agreed guidelines.
- STEP 6** Allow the surface to dry before use.

