


Flower SMART

NITRILE Powder-free Examination Gloves



- ✦ Non sterile, single use gloves
- ✦ Latex free formulation, no natural rubber protein and Powder free, no cornstarch
- ✦ Universal shape suitable for left and right hand (ambidextrous gloves), with beaded cuff
- ✦ Microtextured surface for secure grip
- ✦ Super soft and fatigue free
- ✦ 100 gloves in cardboard box for easy removal of individual gloves
- ✦ Practical box opening to protect the contents from contamination
- ✦ Medical Device Class I in accordance with Regulation (EU) 2017/745
- ✦ PPE category III in accordance with Regulation (EU) 2016/425
- ✦ Hand skin protection against contact with chemicals for a minimum of 30 minutes - Type B
- ✦ Protection against all micro-organisms, including blood-borne viruses
- ✦ Suitable for food contact and use in the food area
- ✦ Colour blue
- ✦ Length min. 240mm
- ✦ AQL ≤ 1.5



Hospital



Diagnostic / Lab



Dental



Patient care



Food contact



Cleaning



Beauty & Wellness

REF	SIZE	PACKING - pcs in box / carton
1323806919	XS ●	100 / 10
1323806920	S ●	100 / 10
1323806921	M ●	100 / 10
1323806922	L ●	100 / 10
1323806923	XL ●	100 / 10

batist@batist.com

+420 491 413311



www.batist.com