

Flower SMART

NITRILE Powder-free Examination Gloves



About

Non-sterile, single-use gloves
 Latex free formulation
 No natural rubber protein and Powder free, no cornstarch

Benefits

Universal shape suitable for left and right hand (ambidextrous gloves), with beaded cuff
 Microtextured surface for secure grip, super soft and fatigue-free
 Protection against all micro-organisms, including blood-borne viruses

Specification

100 gloves in a cardboard box for easy removal of individual gloves
 Practical box opening to protect the contents from contamination
 Colour - blue, Length min. 240mm, AQL ≤ 1.5

Standards Regulations

Medical Device Class I in accordance with Regulation (EU) 2017/745
 PPE category III in accordance with Regulation (EU) 2016/425
 Hand skin protection against contact with chemicals for a min. of 30 minutes - Type B



Hospital



Diagnostic / Lab



Dental



Patient care



Food contact



Cleaning



Beauty & Wellness

| REF | SIZE | PACKING - pcs in box / carton |
|------------|---|-------------------------------|
| 1323806919 | XS ● | 100 / 10 |
| 1323806920 | S ● | 100 / 10 |
| 1323806921 | M ● | 100 / 10 |
| 1323806922 | L ● | 100 / 10 |
| 1323806923 | XL ● | 100 / 10 |

 www.batist.com



E-mail
 CZ info-cz@batist.com
 SK info-sk@batist.com
 PL info-pl@batist.com
 ES batist.espana@batist.com
 PT geral@batist.com
 EX custexport@batist.com



Phone
 +420 491 413311
 +421 244 888 639
 +48 885 520 011
 +34 91 360 97 07
 +351 229 429 520
 +420 491 413 349

batist
 MEDICAL