

immunity

PRODUCT CATALOGUE



**SURGICAL
DRAPES & SETS**

Content

Universal sets.....	4
Surgical sets.....	7
Cardiological sety.....	8
U-split sets.....	10
Orthopaedic sets.....	11
Urology and gynaecology sets.....	15
Ophthalmic sets.....	20
ORL/dental sets.....	22
Other sets.....	24
Speciality drapes.....	25
Utility drapes + hospidrape.....	27
Adhesive drapes.....	28
Aperture drapes and adhesive aperture drapes.....	30
Table covers and components.....	32
Equipment drapes.....	33
Incision film.....	34
Accessories.....	35
Medikit.....	36
Surgical gowns.....	40
Scrub suit set (shirt + pants) and surgical caps.....	41
Surgical gloves latex.....	42
Medical face masks.....	43



DISPOSABLE SURGICAL DRAPE

Disposable surgical drapes

Disposable surgical drapes, gowns, sets and accessories which have been developed to enhance the safety and comfort of the patient, the surgeon and medical staff.


- Safety strong - impermeable
- Extra large absorbent patches
- Very low linting
- Stays well in place, is not slippery
- High puncture resistance - very durable

-  Suitable for procedures with low fluid levels
-  Suitable for procedures with moderate fluid levels
-  Suitable for procedures with high fluid levels

- 1 Large and very efficient absorbent patches, new structure
- 2 Blue/grey nonwoven, very low linting
- 3 A highly elastic multicomponent PE-layer moisture barrier, puncture resistance
- 4 White backsheet with very low linting more comfort to the patient, more strength

-  Absorbent and impermeable fully glued 3-layer material
-  Absorbent and impermeable 2-layer material
-  Absorbent patch - reinforced material
-  Polyethylen foil
-  Elastic material
-  Adhesive edge

Easily identifiable integrated tube holders in drapes with absorbent patch



UNIVERSAL SETS


Universal set ECO

REF 100 003



- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 2 hand towel 30,5 x 34 cm
- 2 adhesive drape 75 x 90 cm (repositionable)
- 1 adhesive drape 180 x 180 cm (repositionable)
- 1 adhesive drape 150 x 240 cm (repositionable)
- 1 instrument table cover 150 x 190 cm

Pack count: 1/6/12




Universal set - Basic

REF 100 005



- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 2 hand towel 30,5 x 34 cm
- 2 adhesive drape 75 x 90 cm
- 1 adhesive drape 180 x 180 cm
- 1 adhesive drape 150 x 240 cm
- 1 instrument table cover 150 x 190 cm

Pack count: 1/6/12




Universal set, reinforcement Medium

REF 100 006



- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 2 adhesive drape 75 x 90 cm, absorbent patches 15 x 50 cm, tube holders
- 1 adhesive drape 180 x 180 cm, absorbent patches 15 x 50 cm, tube holders
- 1 adhesive drape 150 x 240 cm, absorbent patches 15 x 50 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/6/12



UNIVERSAL SETS


Universal set without Mayo, reinforcement Large

REF 100 007



- 1 op-tape 9 x 50 cm
- 4 hand towel 30,5 x 34 cm
- 2 adhesive drape 75 x 90 cm, absorbent patches 25 x 60 cm, tube holders
- 1 adhesive drape 180 x 180 cm, absorbent patches 25 x 60 cm, tube holders
- 1 adhesive drape 150 x 240 cm, absorbent patches 25 x 60 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/6/12




Universal set, reinforcement Medium, 3-layer

REF 100 010



- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 2 adhesive drape 75 x 90 cm, absorbent patches 15 x 50 cm, tube holders
- 1 adhesive drape 180 x 180 cm, absorbent patches 15 x 50 cm, tube holders
- 1 adhesive drape 150 x 240 cm
- 1 instrument table cover 150 x 190 cm

Pack count: 1/5/10




Universal set, reinforcement Large, 3-layer

REF 100 011



- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 2 adhesive drape 75 x 90 cm, absorbent patches 25 x 60 cm, tube holders
- 1 adhesive drape 180 x 180 cm, absorbent patches 25 x 60 cm, tube holders
- 1 adhesive drape 150 x 240 cm, absorbent patches 25 x 60 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/5/10




UNIVERSAL SETS

Universal set pediatric

REF 100 030

- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 2 adhesive drape 75 x 75 cm (repositionable)
- 2 adhesive drape 150 x 150 cm (repositionable)
- 1 instrument table cover 150 x 190 cm

Pack count: 1/7/14




Liposuction set

REF 1230611137



- 1 drape - 2-layers 90 x 120 cm
- 1 drape - 2-layers 150 x 180 cm, adhesive
- 1 drape - 2-layers 100 x 150 cm, adhesive
- 2 drape - 2-layers 75 x 90 cm, adhesive
- 15 lap Sponge prewashed 45 x 45 cm with x-ray tape, loop tied per 5 ks
- 5 gauze ball 30/30 cm
- 1 disposable scalpel blade, No. 11
- 1 disposable scalpel blade, No. 15
- 1 surgical gowns OPEN BACK, size L
- 1 surgical gowns STANDARD, size XL

Pack count: 1/3/6




Laparotomy set 2

REF 1230611140B



- 1 instrument table cover 150 x 190 cm
- 1 laparoscopy drape 230 x 320 cm, with aperture 25 x 30 cm, with Incise film
- 1 op-tape 10 x 50 cm
- 10 lap Sponge prewashed 45 x 45 cm with x-ray tape, loop tied per 5 ks
- 40 non-woven swabs 10 x 20 cm per 20 ks, 4 layers
- 5 gauze ball 30/30 cm, packed in a bag
- 10 gauze ball 7/7 cm, packed in a bag
- 2 hand towel
- 1 instrument table cover 79 x 145 cm absorbent patches
- 1 surgical gowns STANDARD, size L
- 2 surgical gowns STANDARD, size XL

Pack count: 1/5/5



CARDIOLOGICAL SETY


Angiographic set 8

REF 1230611817



- 1 instrument table cover 150 x 190 cm
- 1 angiographic drape 230 x 300 cm, 2 with aperture 7 x 10 cm
- 1 hypodermic syringe 10 ml 2-part Luer
- 1 hypodermic syringe 20 ml 2-part Luer
- 1 needle green 0,8 x 40 mm, 21 G
- 1 disposable scalpel No. 11
- 1 round bowl 500 ml
- 1 round bowl 125 ml
- 1 round bowl 250 ml
- 50 gauze swab 10 x 10 cm, 8 layers
- 1 surgical gowns STANDARD, size XL

Pack count: 1/4/8




Pacemaker implant set 2

REF 1230611821



- 1 drape - 2-layers 75 x 90 cm
- 1 instrument table cover 100 x 150 cm
- 1 drape - 2-layers 150 x 200 cm, adhesive
- 1 drape s with aperture 75 x 100 cm, with aperture 10 cm
- 1 Elastopore+Pad 10 x 15 cm
- 1 bowl 250 ml
- 2 cover for X-ray diameter 90 cm
- 1 surgical gowns SPECIAL, size M
- 1 surgical gowns SPECIAL, size L
- 80 swabs 10 x 10 gauze, 8 layers, folded
- 5 gauze ball 20/19
- 2 syringe 20 ml 2-part Luer
- 2 hand towel
- 1 velcro 2,5 x 30 cm

Pack count: 1/3/6




Basic sterile angiographic set

REF 1230611823



- 1 drape angiographic 230 x 300 cm, 2 with aperture
- 1 instrument table cover 150 x 190 cm
- 4 hand towel

Pack count: 1/10/10




Pacemaker implantation/ablation set sterile

REF 1230611832



- 1 instrument table cover 150 x 190 cm
- 1 drape angiographic 220 x 320 cm, 2 with aperturey Ø 14 cm
- 1 cover for X-ray diameter 140 cm
- 1 surgical gowns STANDARD, size XL
- 2 hand towel

Pack count: 1/8/8



U-SPLIT SETS


U split sheet set M, 3-layer

REF 100 050



- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 adhesive drape 75 x 90 cm
- 1 adhesive split sheet 230 x 260 cm, split 10 x 100 cm, absorbent patches 75 x 140 cm, tube holders
- 1 adhesive drape 150 x 240 cm
- 1 instrument table cover 150 x 190 cm

Pack count: 1/4/8




U split sheet set L, reinforced, 3-layer

REF 100 051



- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 adhesive drape 75 x 90 cm
- 1 adhesive split sheet 230 x 260 cm, split 20 x 100 cm, absorbent patches 75 x 140 cm, tube holders
- 1 adhesive drape 150 x 240 cm, absorbent patches 25 x 60 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/4/8




U split sheet set ECO

REF 100 055



- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 adhesive drape 75 x 90 cm
- 1 adhesive split sheet 230 x 260 cm, split 10 x 85 cm
- 1 adhesive drape 150 x 240 cm
- 1 instrument table cover 150 x 190 cm

Pack count: 1/5/10




Hip set ECO

REF 120 005



- 2 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 1 stockinette 36 x 120 cm
- 1 utility drape 100 x 150 cm
- 1 adhesive drape 150 x 300 cm, absorbent patches 15 x 50 cm, tube holders
- 1 adhesive split sheet 230 x 260 cm, split 20 x 85 cm, absorbent patches 50 x 100 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/4/8




Hip set, reinforcement Large, 3-layer

REF 120 010



- 3 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 elastic stockinette 30 x 120 cm
- 1 elastic support bandage 12 cm x 6 m
- 2 adhesive drape 75 x 90 cm, absorbent patches 25 x 60 cm, tube holders
- 1 adhesive drape 180 x 180 cm, absorbent patches 25 x 60 cm, tube holders
- 1 adhesive drape 150 x 240 cm, absorbent patches 25 x 60 cm, tube holders
- 1 adhesive split sheet 230 x 260 cm, split 20 x 100 cm, absorbent patches 75 x 140 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/3/6




Extremity set

REF 120 030



- 2 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 adhesive drape 75 x 90 cm
- 1 instrument table cover 150 x 190 cm, folded
- 1 extremity drape 230 x 300 cm, elastic aperture 7 cm, absorbent patches 50 x 100 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/4/8



ORTHOPAEDIC SETS


Extremity set

REF 120 031



- 2 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 elastic stockinette 22 x 75 cm
- 1 utility drape 150 x 190 cm
- 1 extremity drape 230 x 300 cm, elastic aperture 7 cm, absorbent patches 50 x 100 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/4/8




Extremity set ECO

REF 120 311



- 2 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 stockinette 36 x 65 cm
- 1 extremity drape 230 x 300 cm, elastic aperture 7 cm, absorbent patches 75 x 80 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/6/12




Hand/foot set, 3-layer

REF 120 032



- 1 Mayo stand cover 79 x 145 cm, reinforced
- 1 extremity drape 230 x 300 cm, elastic aperture 3,5 cm, absorbent patches 50 x 100 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/6/12




Hand/foot small set

REF 1230611363



- 1 instrument table cover 140 x 150 cm
- 1 extremity drape 150 x 180 cm, with aperture 3,5 cm

Pack count: 1/12/24





Hand/foot set 4

REF 120 325

- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 extremity drape 230 x 300 cm, elastic aperture 3,5 cm, absorbent patches 75 x 80 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/6/12




Knee arthroscopy set

REF 120 035

- 2 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 elastic stockinette 22 x 75 cm
- 1 instrument table cover 150 x 190 cm, folded
- 1 knee arthroscopy drape 230 x 320 cm, fluid collection pouch 105 x 90 cm with sieve and suction port, reinforcement wire and 2 elastic apertures Ø 5 cm
- 1 instrument table cover 150 x 190 cm

Pack count: 1/4/8




Knee arthroscopy set ECO

REF 120 355

- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 2 hand towel 30,5 x 34 cm
- 1 stockinette 36 x 65 cm
- 1 extremity drape 230 x 300 cm, elastic aperture 7 cm, absorbent patches 50 x 100 cm, tube holders
- 1 tube cover 15 x 240 cm, perforated tip
- 1 instrument table cover 140 x 150 cm

Pack count: 1/5/10



ORTHOPAEDIC SETS


Knee arthroscopy set

REF 1230611323



- 1 instrument table cover 150 x 190 cm
- 1 extremity drape 3-layer 230 x 300 cm, with aperture 7,5 x 10,5 cm
- 1 tube cover 15 x 240 cm, perforated tip
- 1 op-tape 10 x 50 cm
- 2 hand towel
- 1 stockinette 36 x 65 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced

Pack count: 1/5/10




Shoulder set 2

REF 1230611326



- 1 instrument table cover 150 x 190 cm
- 1 "U" drape 230 x 260 cm, "U" slit 20 x 85 cm
- 2 drape - 2-layers 150 x 180 cm, adhesive
- 1 op-tape 10 x 50 cm
- 2 hand towel
- 1 stockinette 25 x 80 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced

Pack count: 1/4/8




Set a small procedure

REF 1230611352A



- 1 instrument table cover 100 x 150 cm
- 1 drape - 2-layers 100 x 150 cm, adhesive
- 1 drape - 2-layers 90 x 120 cm
- 1 drape - 2-layers 75 x 90 cm, with aperture 6 x 8 cm
- 10 non-woven swabs 10 x 20 cm, 4 layers
- 2 lap Sponge prewashed 45 x 45 cm with RTG
- 5 gauze ball 30/30 cm, packed in a bag
- 2 surgical gowns STANDARD, size XL
- 1 surgical gowns STANDARD, size L

Pack count: 1/4/8




UROLOGY AND GYNAECOLOGY SETS

Gyn/Cyst set

REF 120 050

- 2 leggings 75 x 120 cm
- 1 utility drape 75 x 80 cm
- 1 gyn/cyst drape 75 x 175 cm, aperture 9 x 15 cm
- 1 instrument table cover 150 x 190 cm

Pack count: 1/5/20




Gyn/Cyst set ECO

REF 120 055

- 2 leggings 75 x 120 cm
- 1 gyn/cyst drape 100 x 175 cm, adhesive aperture 9 x 15 cm
- 1 instrument table cover 100 x 150 cm

Pack count: 1/5/20




Gynecology set

REF 120 053

- 3 hand towel 30,5 x 34 cm
- 2 leggings 75 x 120 cm
- 1 drape 75 x 120 cm
- 1 gyn upper drape 115 x 150 cm
- 1 gyn lower drape 140 x 150 cm
- 1 instrument table cover 100 x 150 cm

Pack count: 1/4/16




Gynecology set II

REF 120 051

- 2 hand towel 30,5 x 34 cm
- 1 Mayo stand cover 79 x 145 cm reinforced
- 1 gyn drape 60 x 120 cm with integrated fluid collection pouch 35 x 52 cm with port, adhesive aperture 9 x 12 cm
- 1 gyn drape 270/230 x 260 cm, with aperture 24 x 21 cm with integrated leggings
- 1 velcro tape 2,5 x 30 cm
- 1 instrument table cover 150 x 190 cm

Pack count: 1/5/10




UROLOGY AND GYNAECOLOGY SETS

Vaginal set 1

REF 1230611535

- 1 drape - 2-layers 75 x 90 cm
- 1 instrument table cover 150 x 190 cm
- 1 drape gynecological with integrated sleeves 240 x 275 cm

Pack count: 1/5/10




Vaginal set 2

REF 1230611542B

- 1 instrument table cover 150 x 190 cm
- 1 drape gynecological with integrated sleeves 240 x 275 cm
- 1 op-tape 10 x 50 cm
- 10 non-woven swabs 10 x 20 cm, 4 layers
- 2 lap Sponge prewashed 45 x 45 cm with RTG
- 5 gauze ball 40/40 cm, packed in a bag
- 2 hand towel
- 1 instrument table cover 79 x 145 cm
- 1 surgical gowns STANDARD, size L
- 1 surgical gowns STANDARD, size XL

Pack count: 1/6/6




LAVH 5 set

REF 1230611543A

- 1 instrument table cover 150 x 190 cm
- 1 drape gynecological with integrated sleeves 240/215 x 350 cm
- 1 op-tape 10 x 50 cm
- 1 tube cover 15 x 240 cm, perforated tip
- 2 lap Sponge prewashed 45 x 45 cm with RTG
- 20 non-woven swabs 10 x 20 cm, 4 layers
- 5 lap Sponge prewashed 23 x 33 cm with RTG
- 6 gauze ball 40/40 cm, packed in a bag
- 2 hand towely
- 1 instrument table cover 79 x 145 cm
- 1 surgical gowns STANDARD, size L
- 2 surgical gowns STANDARD, size XL

Pack count: 1/2/4




UROLOGY AND GYNAECOLOGY SETS

Delivery set

REF 120 550

- 6 hand towel 21 x 25 cm
- 1 underpad 60 x 60 cm
- 1 underpad 60 x 90 cm
- 1 baby sheet 56 x 75 cm
- 1 drape (stitching sheet) 75 x 120 cm
- 1 drape (for under the mother) 90 x 150 cm

Pack count: 1/12




C-Section set 4

REF 1230611519

- 1 instrument table cover 150 x 190 cm
- 1 C-section drape 230 x 330 cm
- 1 op-tape 10 x 50 cm
- 1 baby drape 80 x 90 cm
- 2 hand towely
- 1 instrument table cover 79 x 145 cm

Pack count: 1/4/8




C-Section set I

REF 120 060

- 1 op-tape 9 x 50 cm
- 1 návlek na instrumentační stůl zesílený 79 x 145 cm
- 4 hand towel 30,5 x 34 cm
- 1 drape pro dítě 75 x 120 cm
- 1 drape na císařský řez 180/250 x 330 cm, sběrný vak 72 x 72 cm, incizní fólie, with aperture Ø 16 x 18 cm, 2 výpustě
- 1 instrument table cover 150 x 190 cm

Pack count: 1/5/10



UROLOGY AND GYNAECOLOGY SETS


C-Section set II

REF 120 061



- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 baby sheet 75 x 120 cm
- 1 C-section drape 330 x 230 cm, fluid collection pouch 72 x 72 cm, integrated incision film 36 x 36 cm and tube holders, 2 suction ports
- 1 instrument table cover 150 x 190 cm

Pack count: 1/5/10




TUR set

REF 120 072



- 1 op-tape 9 x 50 cm
- 1 velcro tape 2,5 x 30 cm
- 1 hand towel 30,5 x 34 cm
- 1 TUR-drape 175/282 x 172 cm all-in-one, finger cot, apertures 6 cm and 8 cm (adhesive), integrated fluid collection pouch with sieve and suction port
- 1 instrument table cover 150 x 190 cm

Pack count: 1/6/12




Urological set Small

REF 1230611569



- 1 instrument table cover 100 x 150 cm
- 1 drape urological 180 x 300 cm, with aperture 12 cm, adhesive

Pack count: 1/10/20




Abdomino-perineal set ECO

REF 120 485



- 2 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 2 hand towel 30,5 x 34 cm
- 1 Gyn/Abdominal drape 215/240 x 350 cm all-in-one with integrated leg covers and fluid collection pouch, with sieve and suction port, 2 adhesive apertures: top 25 x 30 cm, bottom 9 x 15 cm (with integrated cover drape if apertures not used)
- 1 instrument table cover 150 x 190 cm

Pack count: 1/5/10



OPHTHALMIC SETS



Ophthalmic set

REF 120 080



- 1 hand towel 30,5 x 34 cm
- 1 ophthalmic drape 150 x 140 cm, with fluid collection pouch and integrated incision film with aperture 8 x 10 cm
- 1 instrument table cover 150 x 140 cm

Pack count: 1/7/28





Ophthalmic set II

REF 120 081



- 1 hand towel 30,5 x 34 cm
- 1 ophthalmic drape 150 x 140 cm, with fluid collection pouch and integrated incision film with oval aperture 3,5 x 6,5 cm
- 1 instrument table cover 150 x 140 cm

Pack count: 1/7/28





Ophthalmic set

REF 120 082



- 1 incision film 10 x 12 cm
- 1 ophthalmic drape 150 x 180 cm, 2 fluid collection pouches, oval aperture 5 x 7 cm, 3 reinforcement wires
- 1 instrument table cover 150 x 140 cm

Pack count: 1/7/28




Ophthalmic set outpatient

REF 1230611722



- 1 Utility drape Hospidrape 75 x 90 cm
- 1 drape - 2-layers 50 x 60 cm, with aperture 6 x 8 cm, adhesive

Pack count: 1/70/140




OPHTHALMIC SETS

Ophthalmic set 20

REF 1230611739E

- 1 drape - 2-layers 100 x 150 cm
- 1 drape ophthalmic 120 x 160 cm
- 3 hypodermic syringe 2 ml 2-part Luer
- 1 needle pink 1,2 x 40 18 G
- 1 syringe 2 ml 2-part Luer
- 5 gauze swab 5 x 5 cm, 8 layers
- 1 oval nonwoven ophthalmic compression 7 x 5,5 cm
- 1 round bowl 60 ml
- 2 nonwoven tampon 20/20 cm
- 1 drape - 2-layers 75 x 90 cm
- 1 surgical gowns STANDARD, size M
- 1 surgical gowns STANDARD, size L

Pack count: 1/6/12




Ophthalmic set cataract 2

REF 1230611723

- 1 drape - 2-layers 90 x 120 cm
- 1 drape ophthalmic 120 x 120 cm s with apertureem 7 x 9 cm and 1 holding pocket, incision film
- 2 armrest cover 20 x 55 cm
- 1 sleeve on FACO 56 x 75 cm
- 1 ophthalmic compression 7 x 5,5 cm
- 2 gauze balls 30/30
- 1 needle orange 25 G 0,5 x 25 cm
- 3 syringe 2 ml 2-part Luer
- 1 drape SPUNBOND 80 x 90 cm
- 2 surgical gowns BASIC, size L
- 1 hand towel

Pack count: 1/6/12




ORL/DENTAL SETS

ENT set

REF 120 038

- 2 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 adhesive drape 75 x 90 cm, absorbent patches 15 x 50 cm, tube holders
- 1 utility drape 75 x 80 cm
- 1 adhesive split sheet 150 x 250 cm, split 7 x 70 cm
- 1 instrument table cover 150 x 190 cm

Pack count: 1/5/10




ENT set Small

REF 1230611406

- 1 instrument table cover 90 x 120 cm
- 1 drape 160 x 240 cm, with aperture adhesive 10 cm
- 2 surgical gowns OPEN BACK, size L

Pack count: 1/8/16




ENT set universal Large

REF 1230611416B

- 1 instrument table cover 140 x 150 cm
- 1 drape - 2-layers 180 x 250 cm, with aperture 8 x 10 cm
- 1 lap Sponge prewashed 45 x 45 cm with RTG
- 1 op-tape 10 x 50 cm
- 10 nonwoven compression 10 x 20 cm, 4 layers
- 10 gauze swabs 7,5 x 7,5 cm, 8 layers
- 5 gauze balls 20/19 cm with RTG, packed in a bag
- 1 hand towel
- 1 surgical gowns STANDARD, size L
- 1 surgical gowns STANDARD, size XL

Pack count: 1/10/10




Tympanostomy set

REF 1230611426

- 1 drape - 2-layers 90 x 120 cm
- 1 drape - 2-layers 75 x 90 cm, with aperture 6 x 8 cm
- 1 connecting hose funnel-free end CH30, 2,5 m

Pack count: 1/40/40




ORL-TE set

REF 1230611424

- 1 instrument table cover 100 x 150 cm
- 1 drape dental 100 x 150 cm s with aperture 9,3 x 13 cm
- 10 gauze swab 10 x 10 cm, 8 layers
- 1 syringe 10 ml, 2-part Luer
- 1 scalpel - disposable blade size 15
- 10 gauze balls 20/19 cm, packed in a bag
- 1 suction set Yankauer Ch25, bend 6 mm
- 1 surgical gowns BASIC OPEN BACK, size L

Pack count: 1/8/16




OTHER SETS

Laparoscopy set

REF 120 067

- 2 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 laparoscopy drape 230 x 300 cm, aperture 25 x 20 cm
- 1 fluid collection pouch 30 x 40 cm
- 1 instrument table cover 150 x 190 cm

Pack count: 1/7/14




SPECIALITY DRAPES

Adhesive split sheets

REF 130 010


150 x 250 cm, split 7 x 70 cm, absorbent patch 50 x 50 cm, tube holders
Pack count: 1/8/32



150 x 250 cm
U = 7 x 70 cm

REF 1230612102

230 x 260 cm, split 20 x 85 cm, absorbent patch 50 x 100 cm, tube holders
Pack count: 1/8/16




230 x 260 cm
U = 20 x 85 cm

Hand/foot drape

REF 130 020


150 x 180 cm, elastic aperture Ø 3,5 cm, absorbent patch 50 x 100 cm, tube holders
Pack count: 1/17/34



150 x 180 cm
Ø 3,5 cm

REF 130 024

230 x 300 cm, elastic aperture Ø 3,5 cm, absorbent patch 50 x 100 cm, tube holders
Pack count: 1/8/16




230 x 300 cm
Ø 3,5 cm

Extremity drape

REF 130 025

230 x 300 cm, elastic aperture Ø 7 cm, absorbent patch 50 x 100 cm, tube holders
Pack count: 1/8/16




230 x 300 cm
Ø 7 cm

SPECIALITY DRAPES

Vertical drapes

REF 130 032

240 x 330 cm, incision film 86 x 25 cm, integrated fluid collection pouch with sieve and bottom outlet, 2 tube holders and 2 compartment suction/diathermy pouches, 24 x 24 cm
Pack count: 1/10/20




330 x 240 cm

Laparoscopy drapes

REF 130 066

230 x 300 cm, aperture 20 x 31 cm
Pack count: 1/9/18




230 x 300 cm

Variable aperture drapes

REF 150 012

50 x 75 cm, 2 drapes
Pack count: 1/35/140




50 cm

50 cm

REF 150 013

75 x 100 cm, 2 drapes
Pack count: 1/20/80




75 cm

75 cm

Ophthalmic drapes

REF 1230618103


120 x 120 cm, aperture 7 x 9 cm, with incision foil and fluid collection pouch
Pack count: 1/40/80



120 x 120 cm

REF 1230618104

150 x 150 cm, aperture 7 x 9 cm, with incision foil
Pack count: 1/30/60



150 x 150 cm

Utility DRAPES + HOSPIDRAPE

Utility drapes

REF 140 010

50 x 60 cm

Pack count: 1/45/180

REF 140 015

75 x 90 cm

Pack count: 1/25/100

REF 1230612207

90 x 120 cm

Pack count: 1/70/140

REF 1230612208

100 x 150 cm

Pack count: 1/45/90

REF 1230612209

150 x 180 cm

Pack count: 1/60/60



50 x 60 cm



75 x 90 cm



90 x 120 cm



100 x 150 cm



150 x 180 cm

Hospidrape (Cellulose /PE) - AMBULATORY DRAPES

REF 1230612001

38 x 45 cm

Pack count: 1/90/720



38 x 45 cm

REF 1230612002

50 x 60 cm

Pack count: 1/100/400



50 x 60 cm

REF 1230612003

45 x 75 cm

Pack count: 1/110/440



45 x 75 cm

REF 1230612004

75 x 90 cm

Pack count: 1/120/240



75 x 90 cm

REF 1230612005

90 x 120 cm

Pack count: 1/90/180



90 x 120 cm

ADHESIVE DRAPES

Adhesive drapes - 2-layers

REF 140 055

50 x 60 cm

Pack count: 1/40/160



50 x 60 cm



75 x 90 cm

REF 1230612204


75 x 90 cm

Pack count: 1/100/200


REF 1230612206

150 x 240 cm

Pack count: 1/20/40



150 x 240 cm



90 x 120 cm

REF 1230612212


90 x 120 cm

Pack count: 1/65/130


REF 1230612213

100 x 150 cm

Pack count: 1/40/80



100 x 150 cm



150 x 180 cm

REF 1230612214

150 x 180 cm

Pack count: 1/55/55


ADHESIVE DRAPES

Adhesive drapes - 3-layer

REF 140 060

75 x 90 cm

Pack count: 1/16/64




75 x 90 cm

REF 140 054

75 x 100 cm

Pack count: 1/15/60




75 x 100 cm

REF 140 053

100 x 100 cm

Pack count: 1/10/40




100 x 100 cm

Adhesive drape Medium

REF 140 065

180 x 180 cm

Pack count: 1/8/32




180 x 180 cm

Adhesive drape Large

REF 140 070

150 x 240 cm

Pack count: 1/8/32




150 x 240 cm

APERTURE DRAPES AND ADHESIVE APERTURE DRAPES


Aperture drapes

REF 150 055

50 x 60 cm, Ø 6 x 8 cm
Pack count: 1/45/180



50 x 60 cm
Ø 6 x 8 cm



75 x 90 cm
Ø 6 x 8 cm


REF 150 065

75 x 90 cm, Ø 6 x 8 cm
Pack count: 1/25/100


Adhesive aperture drapes

REF 150 005


50 x 60 cm, Ø 6 x 8 cm
Pack count: 1/45/180



50 x 60 cm
Ø 6 x 8 cm



75 x 80 cm
Ø 6 x 15 cm



75 x 90 cm
Ø 6 x 8 cm

REF 150 010


75 x 80 cm, Ø 6 x 15 cm
Pack count: 1/25/100

REF 150 015


75 x 90 cm, Ø 6 x 8 cm
Pack count: 1/25/100

REF 150 020

120 x 150 cm, Ø 5 x 7 cm
Pack count: 1/20/80



120 x 150 cm
Ø 5 x 7 cm




150 x 180 cm
Ø 5 x 7 cm

REF 150 025


150 x 180 cm, Ø 5 x 7 cm
Pack count: 1/14/56

REF 150 030

150 x 180 cm, Ø 6 x 15 cm
Pack count: 1/12/48



150 x 180 cm
Ø 6 x 15 cm



230 x 250 cm
Ø 12 cm

REF 150 035

230 x 250 cm, Ø 12 cm
Pack count: 1/7/28

APERTURE DRAPES AND ADHESIVE APERTURE DRAPES

REF 1230612310


180 x 280 cm, Ø 14 cm

Pack count: 1/14/28


REF 1230612320

120 x 150 cm, Ø 12 cm

Pack count: 1/40/80



120 x 150 cm
Ø 12 cm



180 x 280 cm
Ø 14 cm


TABLE COVERS AND COMPONENTS

Instrument Table Cover


REF 170 015

100 x 150 cm

Pack count: 1/18/72



100 x 150 cm



150 x 190 cm

REF 170 025

150 x 190 cm


Pack count: 1/15/60

Mayo stand cover reinforced

REF 170 012

79 x 145 cm

Pack count: 1/20/80



79 cm

145 cm

Adhesive OP-tapes

REF 140 095

9 x 50 cm, 1 ks


Pack count: 1/100/400



9 x 50 cm

Tube covers**REF 180 010**


8 x 120 cm, simple tip
Pack count: 1/65/260

**REF 180 015**

15 x 250 cm, simple tip
Pack count: 1/30/120

REF 180 020

15 x 240 cm, perforated tip
Pack count: 1/30/120

**REF 180 025**


17 x 240 cm, elastic tip
Pack count: 1/30/120

Banded Bags & Dome Bags**REF 180 050**

Ø 90 cm
Pack count: 1/18/72

REF 180 055

Ø 140 cm
Pack count: 1/10/40



INCISION FILM

Incizní fólie**REF 1230610101**

15 x 20 cm

Pack count: 1/120/480

REF 1230610102

30 x 20 cm

Pack count: 1/100/200

REF 1230610103

45 x 20 cm

Pack count: 1/15/240

REF 1230610104

40 x 35 cm

Pack count: 1/40/40

REF 1230610105


45 x 50 cm

Pack count: 1/30/30

REF 1230610106

45 x 65 cm


Pack count: 1/30/30



ACCESSORIES

Arm cover**REF 1230613110**

Pack count: 1/100/200

**Leggings****REF 130 090**

75 x 120 cm, 2 pcs/bag

Pack count: 2/18/72

75 x 120 cm
2 ks**Elastic stockinettes****REF 130 095**


30 x 120 cm, size L, + 2 x OP tape 9 x 50 cm

Pack count: 1/10/40

REF 130 096


22 x 75 cm, size S, + 2 x OP tape 9 x 50 cm

Pack count: 1/20/80

**Fluid collection pouches****REF 180 070**

30 x 40 cm, 1 chamber

Pack count: 1/55/220




30 x 40 cm

REF 180 072

30 x 40 cm (30 cm x 15 cm + 25 cm), 2 chambers, suction and diathermy bag

Pack count: 1/55/220




30 x 40 cm

REF 180 073

62 x 50 cm, sieve and outlet

Pack count: 1/11/44




62 x 50 cm

Set 311 for local anesthesia

REF 1230111311

- 1 Utility drape Hospidrape 75 x 90 cm - in this drape the whole set is packed
- 1 aperture drape - 2-layers 75 x 90 cm, aperture 6 x 8 cm, adhesive
- 1 hand towel
- 1 plaster Elastopore+Pad 7 x 5 cm
- 1 syringe 5 ml 2-part Luer
- 1 syringe 10 ml 2-part Luer
- 1 needle black 22 G
- 1 needle pink 18 G
- 1 pean 14 cm
- 1 plastic tweezers green
- 1 scissors Iris 11,5 cm
- 2 medical tray 120 ml
- 5 gauze balls 20/19 cm

Pack count: 1/55/55




Set 160 on a dry wound dressing

REF 1230111160

- 5 gauze swab 7,5 x 7,5 cm, 8 layers
- 1 plastic tweezers green
- 5 gauze balls 20/19 cm

Pack count: 1/80/160




Set 198 to start dialysis

REF 1230111198

- 1 Utility drape Hospidrape 50 x 60 cm
- 2 plaster Softpore 2,5 x 18 cm, 3 strips
- 5 gauze swab 5 x 5 cm, 8 layers
- 5 gauze balls 20/19 cm
- 2 examination gloves, powder free, size M

Pack count: 1/40/80




Set 163 for urinary catheterization

REF 1230111163

- 1 Utility drape Hospidrape 45 x 75 cm
- 1 drape - 2-layers 45 x 75 cm, with aperture 7 cm
- 1 medical tray
- 2 gauze swab 7,5 x 7,5 cm, 8 layers
- 2 plastic tweezers green
- 1 syringe 10 ml 2-part Luer
- 4 gauze balls 20/19 cm
- 2 examination gloves non powdered, latex size M

Pack count: 1/30/30




Set 313 central venous catheterization

REF 1230111313

- 1 Utility drape Hospidrape 75 x 90 cm
- 1 drape - 2-layers 75 x 90 cm, adhesive with aperture 6 x 8 cm
- 1 bowl blue 500 ml
- 2 gauze swab 5 x 5 cm, 8 layers
- 6 nonwoven swabs 7,5 x 7,5 cm, 4 layers
- 1 plaster Dermafoil 10 x 12 cm
- 1 pean 16 cm
- 1 scissors Iris 11,5 cm
- 1 surgical tweezer 12 cm
- 1 clamp Kelly 13 cm
- 1 disposable scalpel No. 11
- 1 syringe 10 ml 2-part Luer
- 1 syringe 20 ml 2-part Luer
- 1 needle black 22 G
- 1 needle pink 18 G
- 1 hand towel
- 5 nonwoven ball 20/20 cm

Pack count: 1/24/24




Set 164 for central venous catheterization

REF 1230111164

- 1 Utility drape Hospidrape 75 x 90 cm
- 1 drape - 2-layers 75 x 90 cm with perforation, adhesive with aperture 6 cm
- 1 plaster Dermafoil 10 x 12 cm
- 1 plastic transparent tray 21 x 14 x 4,5 cm
- 2 gauze swab 5 x 5 cm, 8 layers
- 6 gauze swab 7,5 x 7,5 cm, 8 layers
- 5 gauze balls 30/30 cm
- 1 plastic tweezers green
- 1 scissors Iris 11 cm
- 1 pean 14 cm
- 1 clamp Kelly 13 cm
- 1 disposable scalpel No. 11
- 1 syringe 10 ml 2-part Luer
- 1 syringe 20 ml 2-part Luer
- 1 needle black 22 G
- 1 needle pink 18 G
- 2 examination gloves non powdered, size M

Pack count: 1/20/20




SET 215 for small operations

REF 1230111215

- 1 Utility drape Hospidrape 50 x 60 cm
- 5 lap sponge 45 x 45 cm with X-ray, loop, prewashed
- 20 gauze balls 30/30 cm, packed in a bag
- 5 gauze cuttings 7 x 90 cm

Pack count: 1/30/30




SET 220 for large operations

REF 1230111220

- 1 Utility drape Hospidrape 75 x 90 cm
- 25 lap sponge 45 x 45 cm with X-ray, loop, unwashed
- 20 gauze balls 9/9 cm with X-ray, packed in a bag
- 50 gauze balls 30/30 cm, packed in a bag
- 5 folded gauze cuttings 7 x 25 cm

Pack count: 1/10/10




Set 336 package mammy

REF 1230111336

- 2 rolled surgical cottonwool, 100 cotton, 15 x 15 cm
- 5 cutting gauze folded 17 x 17 cm
- 6 cutting gauze folded 8 x 22 cm
- 2 knitted conforming bandage 20 cm x 10 m
- 4 gauze swab 7,5 x 7,5 cm, 8 layers
- 6 gauze balls 20/19 cm

Pack count: 1/25/25




Set 337 package hand

REF 1230111337

- 3 gauze balls 20/19 cm
- 5 gauze swab 5 x 5 cm, 8 layers
- 10 gauze swab 7,5 x 7,5 cm, 8 layers
- 10 gauze swab 10 x 10 cm, 8 layers
- 1 cotton pad small

Pack count: 1/70/140




Set 536 obstetric

REF 1230111536

- 1 Utility drape Hospidrape 90 x 120 cm
- 1 drape - 2-layers 90 x 120 cm
- 1 baby drape 80 x 90 cm SPUNLACE
- 1 lap sponge 45 x 45 cm with X-ray, prewashed
- 10 gauze balls 40/40

Pack count: 1/23/46




Set 482 for wound dressing

REF 1230111482

- 6 nonwoven ball 20/20 cm
- 1 plastic tweezers green

Pack count: 1/75/150




Set 568 seductive anesthesia for sewing

REF 1230111568

- 1 Utility drape Hospidrape 75 x 90 cm - in this drape the whole set is packed
- 1 aperture drape 75 x 80 cm, adhesive with aperture 6 x 15 cm
- 3 gauze swab 5 x 5 cm, 8 layers
- 5 gauze balls 30/30 cm
- 1 syringe 2 ml 2-part Luer
- 1 syringe 5 ml 2-part Luer
- 1 needle Orange 25 G long
- 2 needle yellow 20 G
- 1 medical tray, 3 parts 227 x 178 x 45 mm
- 1 pean 14 cm

Pack count: 1/30/30




Set 494 drying CHIR

REF 1230111494

- 1 Utility drape Hospidrape 75 x 90 cm
- 1 lap sponge prewashed 45 x 45 cm with X-ray
- 20 gauze balls 30/30 with X-ray, packed in a bag
- 10 gauze balls 20/19 with X-ray, packed in a bag
- 5 gauze balls 9/9 with X-ray, packed in a bag
- 20 gauze swab 10 x 20 cm, 8 layers, folded

Pack count: 1/20/20



SURGICAL GOWNS

SMMMS

SR - BASIC

Open back, for dry operation of short duration



Ref number	Size	Units in carton
250 032	123 cm, L without hand towels	1/20/40
250 033	147 cm, XL without hand towels	1/18/36



SR - Standard

Universal gown, for dry operations



Ref number	Size	Units in carton
250 022	125 cm, M without hand towels	1/20/40
250 023	135 cm, L without hand towels	1/20/40
250 024	145 cm, XL without hand towels	1/18/36
250 025	155 cm, XXL without hand towels	1/18/36
250 027	125 cm, M with 2 hand towels	1/20/40
250 028	135 cm, L with 2 hand towels	1/20/40
250 029	145 cm, XL with 2 hand towels	1/18/36
250 030	155 cm, XXL with 2 hand towels	1/18/36



SR - Special

Reinforced for wet operations


Ref number	Size	Units in carton
250 052	126 cm, M without hand towels	1/20/40
250 053	137 cm, L without hand towels	1/20/40
250 054	147 cm, XL without hand towels	1/18/36
250 055	155 cm, XXL without hand towels	1/18/36
250 057	126 cm, M with 2 hand towels	1/20/40
250 058	137 cm, L with 2 hand towels	1/20/40
250 059	147 cm, XL with 2 hand towels	1/18/36
250 060	155 cm, XXL with 2 hand towels	1/18/36



SURGICAL KITS AND CAPS


Scrub suit set (shirt + pants)

Ref number	Size	Units in carton
1326002661	size S immunity	1/50/50
1326002662	size M immunity	1/50/50
1326002663	size L immunity	1/50/50
1326002664	size XL immunity	1/50/50
1326002665	size XXL immunity	1/50/50



Surgical caps

Ref number	Size	Units in carton
1326002745	Surgical nonwoven cap Bouffant type- white	1/100/1000
1326002746	Surgical nonwoven cap Bouffant type- blue	1/100/1000
1326002747	Surgical nonwoven cap Bouffant type- green	1/100/1000
1326002748	Surgical nonwoven cap Clip type- white	1/100/1000
1326002749	Surgical nonwoven cap Clip type- blue	1/100/1000
1326002750	Surgical nonwoven cap Clip type- green	1/100/1000




SURGICAL GLOVES LATEX

Surgical gloves latex, sterile

Powder-free


Ref number	Size	Units in carton
1323816240	5,5 immunity per 50 pairs	1/50/200
1323816241	6 immunity per 50 pairs	1/50/200
1323816242	6,5 immunity per 50 pairs	1/50/200
1323816243	7 immunity per 50 pairs	1/50/200
1323816244	7,5 immunity per 50 pairs	1/50/200
1323816245	8 immunity per 50 pairs	1/50/200
1323816246	8,5 immunity per 50 pairs	1/50/200
1323816247	9 immunity per 50 pairs	1/50/200



Surgical gloves latex, sterile

Powdered



Ref number	Size	Units in carton
1323816271	6 immunity per 50 pairs	1/50/200
1323816272	6,5 immunity per 50 pairs	1/50/200
1323816273	7 immunity per 50 pairs	1/50/200
1323816274	7,5 immunity per 50 pairs	1/50/200
1323816275	8 immunity per 50 pairs	1/50/200
1323816276	8,5 immunity per 50 pairs	1/50/200
1323816277	9 immunity per 50 pairs	1/50/200




MEDICAL FACE MASKS

NANO FILTER MASK, FACE MASK A OP-MASK


Ref number	Size	Units in carton
1230500529	Nano filter mask B	1/50/1000
1230500519	Nano fiber mask A	1/50/ 750
1230500601	Face mask with ties	1/50/2000
1320100619	OP-MASK A with ties	1/50/2000
1320100629	OP-MASK B with earloops	1/50/2000





immunity

PRODUCT CATALOGUE



BATIST Medical a.s.
Nerudova 309
549 41 Červený Kostelec
Czech Republic

tel.: +420 491 413 311
fax: +420 491 413 394
www.batist.com

