

immunity

PRODUCT CATALOGUE



**SURGICAL
DRAPES & SETS**

Content

Universal sets.....	4
Surgical sets.....	6
Cardiological sety.....	7
U-split sets.....	8
Orthopaedic sets.....	9
Urology and gynaecology sets.....	13
Ophthalmic sets.....	16
ORL/dental sets.....	17
Speciality drapes.....	18
Utility drapes + hospidrape.....	19
Adhesive drapes.....	20
Aperture drapes and adhesive aperture drapes.....	21
Table covers and components.....	22
Equipment drapes.....	23
Incision film.....	24
Accessories.....	25
Medikit.....	26
Surgical gowns and IntimVest.....	29
Surgical caps and surgical gloves latex.....	30
Notes.....	31



DISPOSABLE SURGICAL DRAPE

Disposable surgical drapes

Disposable surgical drapes, gowns, sets and accessories which have been developed to enhance the safety and comfort of the patient, the surgeon and medical staff.


- Safety strong - impermeable
- Extra large absorbent patches
- Very low linting
- Stays well in place, is not slippery
- High puncture resistance - very durable

-  Suitable for procedures with low fluid levels
-  Suitable for procedures with moderate fluid levels
-  Suitable for procedures with high fluid levels

- 1 Large and very efficient absorbent patches, new structure
- 2 Blue/grey nonwoven, very low linting
- 3 A highly elastic multicomponent PE-layer moisture barrier, puncture resistance
- 4 White backsheet with very low linting more comfort to the patient, more strength

-  Absorbent and impermeable fully glued 3-layer material
-  Absorbent and impermeable 2-layer material
-  Absorbent patch - reinforced material
-  Polyethylen foil
-  Elastic material
-  Adhesive edge

Easily identifiable integrated tube holders in drapes with absorbent patch



UNIVERSAL SETS


Universal set - Basic

REF 100 005



- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 2 hand towel 30,5 x 34 cm
- 2 adhesive drape 75 x 90 cm
- 1 adhesive drape 180 x 180 cm
- 1 adhesive drape 150 x 240 cm
- 1 instrument table cover 150 x 190 cm

Pack count: 1/6/12




Universal set, reinforcement Medium

REF 100 006



- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 2 adhesive drape 75 x 90 cm, absorbent patches 15 x 50 cm, tube holders
- 1 adhesive drape 180 x 180 cm, absorbent patches 15 x 50 cm, tube holders
- 1 adhesive drape 150 x 240 cm, absorbent patches 15 x 50 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/6/12



UNIVERSAL SETS


Universal set, reinforcement Medium, 3-layer

REF 100 010



- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 2 adhesive drape 75 x 90 cm, absorbent patches 15 x 50 cm, tube holders
- 1 adhesive drape 180 x 180 cm, absorbent patches 15 x 50 cm, tube holders
- 1 adhesive drape 150 x 240 cm
- 1 instrument table cover 150 x 190 cm

Pack count: 1/5/10




Universal set, reinforcement Large, 3-layer

REF 100 011



- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 2 adhesive drape 75 x 90 cm, absorbent patches 25 x 60 cm, tube holders
- 1 adhesive drape 180 x 180 cm, absorbent patches 25 x 60 cm, tube holders
- 1 adhesive drape 150 x 240 cm, absorbent patches 25 x 60 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/5/10



SURGICAL SETS


Liposuction set

REF 1230611137



- 1 drape - 2-layers 90 x 120 cm
- 1 drape - 2-layers 150 x 180 cm, adhesive
- 1 drape - 2-layers 100 x 150 cm, adhesive
- 2 drape - 2-layers 75 x 90 cm, adhesive
- 15 lap Sponge prewashed 45 x 45 cm with x-ray tape, loop tied per 5 ks
- 5 gauze ball 30/30 cm
- 1 disposable scalpel blade, No. 11
- 1 disposable scalpel blade, No. 15
- 1 surgical gowns OPEN BACK, size L
- 1 surgical gowns STANDARD, size XL

Pack count: 1/3/6




Laparotomy set 2

REF 1230611140B



- 1 instrument table cover 150 x 190 cm
- 1 laparoscopy drape 230 x 320 cm, with aperture 25 x 30 cm, with Incise film
- 1 op-tape 10 x 50 cm
- 10 lap Sponge prewashed 45 x 45 cm with x-ray tape, loop tied per 5 ks
- 40 non-woven swabs 10 x 20 cm per 20 ks, 4 layers
- 5 gauze ball 30/30 cm, packed in a bag
- 10 gauze ball 7/7 cm, packed in a bag
- 2 hand towel
- 1 instrument table cover 79 x 145 cm absorbent patches
- 1 surgical gowns STANDARD, size L
- 2 surgical gowns STANDARD, size XL

Pack count: 1/5/5



CARDIOLOGICAL SET


Pacemaker implant set 2

REF 1230611821



- 1 drape - 2-layers 75 x 90 cm
- 1 instrument table cover 100 x 150 cm
- 1 drape - 2-layers 150 x 200 cm, adhesive
- 1 drape s with apertureem 75 x 100 cm, with aperture 10 cm
- 1 Elastpore+Pad 10 x 15 cm
- 1 bowl 250 ml
- 2 cover for X-ray diameter 90 cm
- 1 surgical gowns SPECIAL, size M
- 1 surgical gowns SPECIAL, size L
- 80 swabs 10 x 10 gauze, 8 layers, folded
- 5 gauze ball 20/19
- 2 syringe 20 ml 2-part Luer
- 2 hand towel
- 1 velcro 2,5 x 30 cm

Pack count: 1/3/6




Pacemaker implantation/ablation set sterile

REF 1230611832



- 1 instrument table cover 150 x 190 cm
- 1 drape angiographic 220 x 320 cm, 2 with aperturey ø 14 cm
- 1 cover for X-ray diameter 140 cm
- 1 surgical gowns STANDARD, size XL
- 2 hand towel

Pack count: 1/8/8



U-SPLIT SETS


U split sheet set M, 3-layer

REF 100 050



- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 adhesive drape 75 x 90 cm
- 1 adhesive split sheet 230 x 260 cm, split 10 x 100 cm, absorbent patches 75 x 140 cm, tube holders
- 1 adhesive drape 150 x 240 cm
- 1 instrument table cover 150 x 190 cm

Pack count: 1/4/8




U split sheet set L, reinforced, 3-layer

REF 100 051



- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 adhesive drape 75 x 90 cm
- 1 adhesive split sheet 230 x 260 cm, split 20 x 100 cm, absorbent patches 75 x 140 cm, tube holders
- 1 adhesive drape 150 x 240 cm, absorbent patches 25 x 60 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/4/8




U split sheet set ECO

REF 100 055



- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 adhesive drape 75 x 90 cm
- 1 adhesive split sheet 230 x 260 cm, split 10 x 85 cm
- 1 adhesive drape 150 x 240 cm
- 1 instrument table cover 150 x 190 cm

Pack count: 1/5/10




Hip set ECO

REF 120 005



- 2 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 1 stockinette 36 x 120 cm
- 1 utility drape 100 x 150 cm
- 1 adhesive drape 150 x 300 cm, absorbent patches 15 x 50 cm, tube holders
- 1 adhesive split sheet 230 x 260 cm, split 20 x 85 cm, absorbent patches 50 x 100 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/4/8




Hip set, reinforcement Large, 3-layer

REF 120 010



- 3 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 elastic stockinette 30 x 120 cm
- 1 elastic support bandage 12 cm x 6 m
- 2 adhesive drape 75 x 90 cm, absorbent patches 25 x 60 cm, tube holders
- 1 adhesive drape 180 x 180 cm, absorbent patches 25 x 60 cm, tube holders
- 1 adhesive drape 150 x 240 cm, absorbent patches 25 x 60 cm, tube holders
- 1 adhesive split sheet 230 x 260 cm, split 20 x 100 cm, absorbent patches 75 x 140 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/3/6




Extremity set

REF 120 030



- 2 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 adhesive drape 75 x 90 cm
- 1 instrument table cover 150 x 190 cm, folded
- 1 extremity drape 230 x 300 cm, elastic aperture 7 cm, absorbent patches 50 x 100 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/4/8



ORTHOPAEDIC SETS


Extremity set

REF 120 031



- 2 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 elastic stockinette 22 x 75 cm
- 1 utility drape 150 x 190 cm
- 1 extremity drape 230 x 300 cm, elastic aperture 7 cm, absorbent patches 50 x 100 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/4/8




Extremity set ECO

REF 120 311



- 2 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 stockinette 36 x 65 cm
- 1 extremity drape 230 x 300 cm, elastic aperture 7 cm, absorbent patches 75 x 80 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/6/12




Hand/foot set, 3-layer

REF 120 032



- 1 Mayo stand cover 79 x 145 cm, reinforced
- 1 extremity drape 230 x 300 cm, elastic aperture 3,5 cm, absorbent patches 50 x 100 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/6/12





Hand/foot set 4

REF 120 325

- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 extremity drape 230 x 300 cm, elastic aperture 3,5 cm, absorbent patches 75 x 80 cm, tube holders
- 1 instrument table cover 150 x 190 cm

Pack count: 1/6/12




Knee arthroscopy set

REF 120 035

- 2 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 elastic stockinette 22 x 75 cm
- 1 instrument table cover 150 x 190 cm, folded
- 1 knee arthroscopy drape 230 x 320 cm, fluid collection pouch 105 x 90 cm with sieve and suction port, reinforcement wire and 2 elastic apertures Ø 5 cm
- 1 instrument table cover 150 x 190 cm

Pack count: 1/4/8




Knee arthroscopy set ECO

REF 120 355

- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 2 hand towel 30,5 x 34 cm
- 1 stockinette 36 x 65 cm
- 1 extremity drape 230 x 300 cm, elastic aperture 7 cm, absorbent patches 50 x 100 cm, tube holders
- 1 tube cover 15 x 240 cm, perforated tip
- 1 instrument table cover 140 x 150 cm

Pack count: 1/5/10



ORTHOPAEDIC SETS


Knee arthroscopy set

REF 1230611323



- 1 instrument table cover 150 x 190 cm
- 1 extremity drape 3-layer 230 x 300 cm, with aperture 7,5 x 10,5 cm
- 1 tube cover 15 x 240 cm, perforated tip
- 1 op-tape 10 x 50 cm
- 2 hand towel
- 1 stockinette 36 x 65 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced

Pack count: 1/5/10




Set a small procedure

REF 1230611352A



- 1 instrument table cover 100 x 150 cm
- 1 drape - 2-layers 100 x 150 cm, adhesive
- 1 drape - 2-layers 90 x 120 cm
- 1 drape - 2-layers 75 x 90 cm, with aperture 6 x 8 cm
- 10 non-woven swabs 10 x 20 cm, 4 layers
- 2 lap Sponge prewashed 45 x 45 cm with RTG
- 5 gauze ball 30/30 cm, packed in a bag
- 2 surgical gowns STANDARD, size XL
- 1 surgical gowns STANDARD, size L

Pack count: 1/4/8




UROLOGY AND GYNAECOLOGY SETS

Gyn/Cyst set ECO

REF 120 055

- 2 leggings 75 x 120 cm
- 1 gyn/cyst drape 100 x 175 cm, adhesive aperture 9 x 15 cm
- 1 instrument table cover 100 x 150 cm

Pack count: 1/5/20




Vaginal set 1

REF 1230611535

- 1 drape - 2-layers 75 x 90 cm
- 1 instrument table cover 150 x 190 cm
- 1 drape gynecological with integrated sleeves 240 x 275 cm

Pack count: 1/5/10




Vaginal set 2

REF 1230611542B

- 1 instrument table cover 150 x 190 cm
- 1 drape gynecological with integrated sleeves 240 x 275 cm
- 1 op-tape 10 x 50 cm
- 10 non-woven swabs 10 x 20 cm, 4 layers
- 2 lap Sponge prewashed 45 x 45 cm with RTG
- 5 gauze ball 40/40 cm, packed in a bag
- 2 hand towel
- 1 instrument table cover 79 x 145 cm
- 1 surgical gowns STANDARD, size L
- 1 surgical gowns STANDARD, size XL

Pack count: 1/6/6




Delivery set

REF 120 550

- 6 hand towel 21 x 25 cm
- 1 underpad 60 x 60 cm
- 1 underpad 60 x 90 cm
- 1 baby sheet 56 x 75 cm
- 1 drape (stitching sheet) 75 x 120 cm
- 1 drape (for under the mother) 90 x 150 cm

Pack count: 1/12




UROLOGY AND GYNAECOLOGY SETS

C-Section set 4

REF 1230611519

- 1 instrument table cover 150 x 190 cm
- 1 C-section drape 230 x 330 cm
- 1 op-tape 10 x 50 cm
- 1 baby drape 80 x 90 cm
- 2 hand towels
- 1 instrument table cover 79 x 145 cm

Pack count: 1/4/8




C-Section set I

REF 120 060

- 1 op-tape 9 x 50 cm
- 1 návlek na instrumentační stůl zesílený 79 x 145 cm
- 4 hand towel 30,5 x 34 cm
- 1 drape pro dítě 75 x 120 cm
- 1 drape na císařský řez 180/250 x 330 cm, sběrný vak 72 x 72 cm, incizní fólie, with aperture Ø 16 x 18 cm, 2 výpustě
- 1 instrument table cover 150 x 190 cm

Pack count: 1/5/10




C-Section set II

REF 120 061

- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 4 hand towel 30,5 x 34 cm
- 1 baby sheet 75 x 120 cm
- 1 C-section drape 330 x 230 cm, fluid collection pouch 72 x 72 cm, integrated incision film 36 x 36 cm and tube holders, 2 suction ports
- 1 instrument table cover 150 x 190 cm

Pack count: 1/5/10





UROLOGY AND GYNAECOLOGY SETS

TUR set

REF 120 072

- 1 op-tape 9 x 50 cm
- 1 velcro tape 2,5 x 30 cm
- 1 hand towel 30,5 x 34 cm
- 1 TUR-drape 175/282 x 172 cm all-in-one, finger cot, apertures 6 cm and 8 cm (adhesive), integrated fluid collection pouch with sieve and suction port
- 1 instrument table cover 150 x 190 cm

Pack count: 1/6/12




Urological set Small

REF 1230611569

- 1 instrument table cover 100 x 150 cm
- 1 drape urological 180 x 300 cm, with aperture 12 cm, adhesive

Pack count: 1/10/20




Abdomino-perineal set ECO

REF 120 485

- 2 op-tape 9 x 50 cm
- 1 Mayo stand cover 79 x 145 cm, reinforced
- 2 hand towel 30,5 x 34 cm
- 1 Gyn/Abdominal drape 215/240 x 350 cm all-in-one with integrated leg covers and fluid collection pouch, with sieve and suction port, 2 adhesive apertures: top 25 x 30 cm, bottom 9 x 15 cm (with integrated cover drape if apertures not used)
- 1 instrument table cover 150 x 190 cm

Pack count: 1/5/10



OPHTHALMIC SETS



Ophthalmic set

REF 120 080



- 1 hand towel 30,5 x 34 cm
- 1 ophthalmic drape 150 x 140 cm, with fluid collection pouch and integrated incision film with aperture 8 x 10 cm
- 1 instrument table cover 150 x 140 cm

Pack count: 1/7/28





Ophthalmic set

REF 120 082



- 1 incision film 10 x 12 cm
- 1 ophthalmic drape 150 x 180 cm, 2 fluid collection pouches, oval aperture 5 x 7 cm, 3 reinforcement wires
- 1 instrument table cover 150 x 140 cm

Pack count: 1/7/28




Ophthalmic set outpatient

REF 1230611722



- 1 Utility drape Hospidrape 75 x 90 cm
- 1 drape - 2-layers 50 x 60 cm, with aperture 6 x 8 cm, adhesive

Pack count: 1/70/140




Ophthalmic set cataract 2

REF 1230611723



- 1 drape - 2-layers 90 x 120 cm
- 1 drape ophthalmic 120 x 120 cm s with aperture 7 x 9 cm and 1 holding pocket, incision film
- 2 armrest cover 20 x 55 cm
- 1 sleeve on FACO 56 x 75 cm
- 1 ophthalmic compression 7 x 5,5 cm
- 2 gauze balls 30/30
- 1 needle orange 25 G 0,5 x 25 cm
- 3 syringe 2 ml 2-part Luer
- 1 drape SPUNBOND 80 x 90 cm
- 2 surgical gowns BASIC, size L
- 1 hand towel

Pack count: 1/6/12




ORL/DENTAL SETS

ENT set Small

REF 1230611406

- 1 instrument table cover 90 x 120 cm
- 1 drape 160 x 240 cm, with aperture adhesive 10 cm
- 2 surgical gowns OPEN BACK, size L

Pack count: 1/8/16




ENT set universal Large

REF 1230611416B

- 1 instrument table cover 140 x 150 cm
- 1 drape - 2-layers 180 x 250 cm, with aperture 8 x 10 cm
- 1 lap Sponge prewashed 45 x 45 cm with RTG
- 1 op-tape 10 x 50 cm
- 10 nonwoven compression 10 x 20 cm, 4 layers
- 10 gauze swabs 7,5 x 7,5 cm, 8 layers
- 5 gauze balls 20/19 cm with RTG, packed in a bag
- 1 hand towel
- 1 surgical gowns STANDARD, size L
- 1 surgical gowns STANDARD, size XL

Pack count: 1/10/10




Tympanostomy set

REF 1230611426

- 1 drape - 2-layers 90 x 120 cm
- 1 drape - 2-layers 75 x 90 cm, with aperture 6 x 8 cm
- 1 connecting hose funnel-free end CH30, 2,5 m

Pack count: 1/40/40




ORL-TE set

REF 1230611424

- 1 instrument table cover 100 x 150 cm
- 1 drape dental 100 x 150 cm s with apertureem 9,3 x 13 cm
- 10 gauze swab 10 x 10 cm, 8 layers
- 1 syringe 10 ml, 2-part Luer
- 1 scalpel - disposable blade size 15
- 10 gauze balls 20/19 cm, packed in a bag
- 1 suction set Yankauer Ch25, bend 6 mm
- 1 surgical gowns BASIC OPEN BACK, size L

Pack count: 1/8/16




SPECIALITY DRAPES

Adhesive split sheets

REF 130 010

150 x 250 cm, split 7 x 70 cm, absorbent patch 50 x 50 cm, tube holders
Pack count: 1/8/32




150 x 250 cm
U = 7 x 70 cm

Hand/foot drape

REF 130 020

150 x 180 cm, elastic aperture Ø 3,5 cm, absorbent patch 50 x 100 cm,
tube holders
Pack count: 1/17/34




150 x 180 cm
Ø 3,5 cm

Ophthalmic drapes

REF 1230618103

120 x 120 cm, aperture 7 x 9 cm, with incision foil
and fluid collection pouch
Pack count: 1/40/80



120 x 120 cm

Utility DRAPES + HOSPIDRAPE

Utility drapes

REF 140 010

50 x 60 cm

Pack count: 1/45/180

REF 140 015

75 x 90 cm

Pack count: 1/25/100

REF 1230612207

90 x 120 cm

Pack count: 1/70/140

REF 1230612208

100 x 150 cm

Pack count: 1/45/90

REF 1230612209

150 x 180 cm

Pack count: 1/60/60



50 x 60 cm



75 x 90 cm



90 x 120 cm



100 x 150 cm



150 x 180 cm

Hospidrape (Cellulose /PE) - AMBULATORY DRAPES

REF 1230612001

38 x 45 cm

Pack count: 1/90/720



38 x 45 cm

REF 1230612002

50 x 60 cm

Pack count: 1/100/400



50 x 60 cm

REF 1230612003

45 x 75 cm

Pack count: 1/110/440



45 x 75 cm

REF 1230612004

75 x 90 cm

Pack count: 1/120/240



75 x 90 cm

REF 1230612005

90 x 120 cm

Pack count: 1/90/180



90 x 120 cm

ADHESIVE DRAPES

Adhesive drapes - 2-layers

REF 1230612204

75 x 90 cm

Pack count: 1/100/200




75 x 90 cm


REF 1230612206

150 x 240 cm

Pack count: 1/20/40



150 x 240 cm




90 x 120 cm


REF 1230612212

90 x 120 cm

Pack count: 1/65/130



100 x 150 cm



150 x 180 cm

REF 1230612213

100 x 150 cm

Pack count: 1/40/80

REF 1230612214

150 x 180 cm


Pack count: 1/55/55

Adhesive drapes - 3-layer

REF 140 060

75 x 90 cm

Pack count: 1/16/64




75 x 90 cm

APERTURE DRAPES AND ADHESIVE APERTURE DRAPES


Aperture drapes

REF 150 055

50 x 60 cm, Ø 6 x 8 cm
Pack count: 1/45/180



50 x 60 cm
Ø 6 x 8 cm



75 x 90 cm
Ø 6 x 8 cm


REF 150 065

75 x 90 cm, Ø 6 x 8 cm
Pack count: 1/25/100


Adhesive aperture drapes

REF 150 005


50 x 60 cm, Ø 6 x 8 cm
Pack count: 1/45/180



50 x 60 cm
Ø 6 x 8 cm



75 x 80 cm
Ø 6 x 15 cm



75 x 90 cm
Ø 6 x 8 cm

REF 150 010


75 x 80 cm, Ø 6 x 15 cm
Pack count: 1/25/100

REF 150 015


75 x 90 cm, Ø 6 x 8 cm
Pack count: 1/25/100

REF 150 020

120 x 150 cm, Ø 5 x 7 cm
Pack count: 1/20/80




120 x 150 cm
Ø 5 x 7 cm



150 x 180 cm
Ø 6 x 15 cm

REF 150 030

150 x 180 cm, Ø 6 x 15 cm
Pack count: 1/12/48



120 x 150 cm
Ø 12 cm

REF 1230612320

120 x 150 cm, Ø 12 cm
Pack count: 1/40/80


TABLE COVERS AND COMPONENTS

Instrument Table Cover


REF 170 015

100 x 150 cm

Pack count: 1/18/72



100 x 150 cm



150 x 190 cm

REF 170 025

150 x 190 cm


Pack count: 1/15/60

Mayo stand cover reinforced

REF 170 012

79 x 145 cm

Pack count: 1/20/80



79 cm

145 cm

Adhesive OP-tapes

REF 140 095

9 x 50 cm, 1 ks


Pack count: 1/100/400



9 x 50 cm

Tube covers**REF 180 010**


8 x 120 cm, simple tip
Pack count: 1/65/260

**REF 180 015**

15 x 250 cm, simple tip
Pack count: 1/30/120

REF 180 020

15 x 240 cm, perforated tip
Pack count: 1/30/120

**REF 180 025**


17 x 240 cm, elastic tip
Pack count: 1/30/120

Banded Bags & Dome Bags**REF 180 050**

Ø 90 cm
Pack count: 1/18/72

REF 180 055

Ø 140 cm
Pack count: 1/10/40



INCISION FILM

Incizní fólie**REF 1230610102**


30 x 20 cm

Pack count: 1/100/200

REF 1230610103

45 x 20 cm

Pack count: 1/15/240

**REF 1230610105**

45 x 50 cm

Pack count: 1/30/30

REF 1230610106


45 x 65 cm

Pack count: 1/30/30

ACCESSORIES


Arm cover**REF 1230613110**

Pack count: 1/100/200

**Elastic stockinettes****REF 130 096**


22 x 75 cm, size S, + 2 x OP tape 9 x 50 cm

Pack count: 1/20/80

**Fluid collection pouches****REF 180 070**


30 x 40 cm, 1 chamber

Pack count: 1/55/220

**REF 180 072**


30 x 40 cm (30 cm x 15 cm + 25 cm), 2 chambers, suction and diathermy bag

Pack count: 1/55/220

**REF 180 073**

62 x 50 cm, sieve and outlet

Pack count: 1/11/44




Set 311 for local anesthesia

REF 1230111311

- 1 Utility drape Hospidrape 75 x 90 cm - in this drape the whole set is packed
- 1 aperture drape - 2-layers 75 x 90 cm, aperture 6 x 8 cm, adhesive
- 1 hand towel
- 1 plaster Elastopore+Pad 7 x 5 cm
- 1 syringe 5 ml 2-part Luer
- 1 syringe 10 ml 2-part Luer
- 1 needle black 22 G
- 1 needle pink 18 G
- 1 pean 14 cm
- 1 plastic tweezers green
- 1 scissors Iris 11,5 cm
- 2 medical tray 120 ml
- 5 gauze balls 20/19 cm

Pack count: 1/55/55




Set 160 on a dry wound dressing

REF 1230111160

- 5 gauze swab 7,5 x 7,5 cm, 8 layers
- 1 plastic tweezers green
- 5 gauze balls 20/19 cm

Pack count: 1/80/160




Set 198 to start dialysis

REF 1230111198

- 1 Utility drape Hospidrape 50 x 60 cm
- 2 plaster Softpore 2,5 x 18 cm, 3 strips
- 5 gauze swab 5 x 5 cm, 8 layers
- 5 gauze balls 20/19 cm
- 2 examination gloves, powder free, size M

Pack count: 1/40/80




Set 163 for urinary catheterization

REF 1230111163

- 1 Utility drape Hospidrape 45 x 75 cm
- 1 drape - 2-layers 45 x 75 cm, with aperture 7 cm
- 1 medical tray
- 2 gauze swab 7,5 x 7,5 cm, 8 layers
- 2 plastic tweezers green
- 1 syringe 10 ml 2-part Luer
- 4 gauze balls 20/19 cm
- 2 examination gloves non powdered, latex size M

Pack count: 1/30/30




Set 313 central venous catheterization

REF 1230111313

- 1 Utility drape Hospidrape 75 x 90 cm
- 1 drape - 2-layers 75 x 90 cm, adhesive with aperture 6 x 8 cm
- 1 bowl blue 500 ml
- 2 gauze swab 5 x 5 cm, 8 layers
- 6 nonwoven swabs 7,5 x 7,5 cm, 4 layers
- 1 plaster Dermafoil 10 x 12 cm
- 1 pean 16 cm
- 1 scissors Iris 11,5 cm
- 1 surgical tweezer 12 cm
- 1 clamp Kelly 13 cm
- 1 disposable scalpel No. 11
- 1 syringe 10 ml 2-part Luer
- 1 syringe 20 ml 2-part Luer
- 1 needle black 22 G
- 1 needle pink 18 G
- 1 hand towel
- 5 nonwoven ball 20/20 cm

Pack count: 1/24/24




Set 164 for central venous catheterization

REF 1230111164

- 1 Utility drape Hospidrape 75 x 90 cm
- 1 drape - 2-layers 75 x 90 cm with perforation, adhesive with aperture 6 cm
- 1 plaster Dermafoil 10 x 12 cm
- 1 plastic transparent tray 21 x 14 x 4,5 cm
- 2 gauze swab 5 x 5 cm, 8 layers
- 6 gauze swab 7,5 x 7,5 cm, 8 layers
- 5 gauze balls 30/30 cm
- 1 plastic tweezers green
- 1 scissors Iris 11 cm
- 1 pean 14 cm
- 1 clamp Kelly 13 cm
- 1 disposable scalpel No. 11
- 1 syringe 10 ml 2-part Luer
- 1 syringe 20 ml 2-part Luer
- 1 needle black 22 G
- 1 needle pink 18 G
- 2 examination gloves non powdered, size M

Pack count: 1/20/20




SET 215 for small operations

REF 1230111215

- 1 Utility drape Hospidrape 50 x 60 cm
- 5 lap sponge 45 x 45 cm with X-ray, loop, prewashed
- 20 gauze balls 30/30 cm, packed in a bag
- 5 gauze cuttings 7 x 90 cm

Pack count: 1/30/30




Set 336 package mammy

REF 1230111336

- 2 rolled surgical cottonwool, 100 cotton, 15 x 15 cm
- 5 cutting gauze folded 17 x 17 cm
- 6 cutting gauze folded 8 x 22 cm
- 2 knitted conforming bandage 20 cm x 10 m
- 4 gauze swab 7,5 x 7,5 cm, 8 layers
- 6 gauze balls 20/19 cm

Pack count: 1/25/25




Set 337 package hand

REF 1230111337

- 3 gauze balls 20/19 cm
- 5 gauze swab 5 x 5 cm, 8 layers
- 10 gauze swab 7,5 x 7,5 cm, 8 layers
- 10 gauze swab 10 x 10 cm, 8 layers
- 1 cotton pad small

Pack count: 1/70/140




Set 482 for wound dressing

REF 1230111482

- 6 nonwovwn ball 20/20 cm
- 1 plastic tweezers green

Pack count: 1/75/150




Set 494 drying CHIR

REF 1230111494

- 1 Utility drape Hospidrape 75 x 90 cm
- 1 lap sponge prewashed 45 x 45 cm with X-ray
- 20 gauze balls 30/30 with X-ray, packed in a bag
- 10 gauze balls 20/19 with X-ray, packed in a bag
- 5 gauze balls 9/9 with X-ray, packed in a bag
- 20 gauze swab 10 x 20 cm, 8 layers, folded

Pack count: 1/20/20





SURGICAL GOWNS SMMSS AND INTIMVEST

SR - Standard

Universal gown, for dry operations



Ref number	Size	Units in carton
250 022B	125 cm, M without hand towels	1/20/40
250 023B	135 cm, L without hand towels	1/20/40
250 024B	145 cm, XL without hand towels	1/18/36
250 025B	155 cm, XXL without hand towels	1/18/36
250 027B	125 cm, M with 2 hand towels	1/20/40
250 028B	135 cm, L with 2 hand towels	1/20/40
250 029B	145 cm, XL with 2 hand towels	1/18/36
250 030B	155 cm, XXL with 2 hand towels	1/18/36



SR - Special


Reinforced for wet operations

Ref number	Size	Units in carton
250 052B	126 cm, M without hand towels	1/20/40
250 053B	137 cm, L without hand towels	1/20/40
250 054B	147 cm, XL without hand towels	1/18/36
250 055B	155 cm, XXL without hand towels	1/18/36
250 057B	126 cm, M with 2 hand towels	1/20/40
250 058B	137 cm, L with 2 hand towels	1/20/40
250 059B	147 cm, XL with 2 hand towels	1/18/36
250 060	155 cm, XXL with 2 hand towels	1/18/36



IntimVest



Ref number	Size	Units in carton
1230500313	universal	6 x 15 pcs



SURGICAL CAPS AND SURGICAL GLOVES LATEX

Surgical caps


Ref number		Units in carton
1326002746	Surgical nonwoven cap Bouffant type- blue	1/100/1000
1326002749	Surgical nonwoven cap Clip type- blue	1/100/1000



Surgical gloves latex, sterile

Powder-free


Ref number	Size	Units in carton
1323816240	5,5 immunity per 50 pairs	1/50/200
1323816241	6 immunity per 50 pairs	1/50/200
1323816242	6,5 immunity per 50 pairs	1/50/200
1323816243	7 immunity per 50 pairs	1/50/200
1323816244	7,5 immunity per 50 pairs	1/50/200
1323816245	8 immunity per 50 pairs	1/50/200
1323816246	8,5 immunity per 50 pairs	1/50/200
1323816247	9 immunity per 50 pairs	1/50/200



Surgical gloves latex, sterile


Powdered

Ref number	Size	Units in carton
1323816271	6 immunity per 50 pairs	1/50/200
1323816272	6,5 immunity per 50 pairs	1/50/200
1323816273	7 immunity per 50 pairs	1/50/200
1323816274	7,5 immunity per 50 pairs	1/50/200
1323816275	8 immunity per 50 pairs	1/50/200
1323816276	8,5 immunity per 50 pairs	1/50/200
1323816277	9 immunity per 50 pairs	1/50/200




Notes

[illegible]



immunity

PRODUCT CATALOGUE



BATIST Medical a.s.
Nerudova 309
549 41 Červený Kostelec
Czech Republic

tel.: +420 491 413 311
fax: +420 491 413 394
www.batist.com

