



immunity

PRODUCT CATALOGUE



**SURGICAL
DRAPES &
SETS**



DISPOSABLE SURGICAL DRAPES

Disposable surgical drapes

Disposable surgical drapes, gowns, sets and accessories which have been developed to enhance the safety and comfort of the patient, the surgeon and medical staff.

- Safety strong - impermeable
- Extra large absorbent patches
- Very low linting
- Stays well in place, is not slippery
- High puncture resistance - very durable



Suitable for procedures with low fluid levels



Suitable for procedures with moderate fluid levels



Suitable for procedures with high fluid levels

- 1 Large and very efficient absorbent patches, new structure
- 2 Blue/grey nonwoven, very low linting
- 3 A highly elastic multicomponent PE-layer moisture barrier, puncture resistance
- 4 White backsheet with very low linting more comfort to the patient, more strength



Absorbent and impermeable fully glued 3-layer material



Absorbent and impermeable 2-layer material



Absorbent patch - reinforced material



Polyethylen foil




Elastic material



Adhesive edge

Easily identifiable integrated tube holders in drapes with absorbent patch



Complies to EN 13795


Universal set - Basic



REF 100 005 / 1320611103

Number of packages:
1/6/12

quantity	size	description
1	9 x 50 cm	op-tape
1	79 x 145 cm	Mayo stand cover, reinforced
2	30,5 x 34 cm	hand towel
2	75 x 90 cm	adhesive drape
1	180 x 180 cm	adhesive drape
1	150 x 240 cm	adhesive drape
1	150 x 190 cm	instrument table cover




Universal set, reinforcement



REF 100 006 / 1320611101

Number of packages:
1/6/12

quantity	size	description
1	9 x 50 cm	op-tape
1	79 x 145 cm	Mayo stand cover, reinforced
4	30,5 x 34 cm	hand towel
2	75 x 90 cm	adhesive drape, absorbent patches 15 x 50 cm, tube holders
1	180 x 180 cm	adhesive drape, absorbent patches 15 x 50 cm, tube holders
1	150 x 240 cm	adhesive drape, absorbent patches 15 x 50 cm, tube holders
1	150 x 190 cm	instrument table cover




U split sheet set ECO



REF 100 055 / 1320611106

Number of packages:
1/5/10


quantity	size	description
1	9 x 50 cm	op-tape
1	79 x 145 cm	Mayo stand cover, reinforced
4	30,5 x 34 cm	hand towel
1	75 x 90 cm	adhesive drape
1	230 x 260 cm	adhesive split sheet, split 10 x 85 cm
1	150 x 240 cm	adhesive drape
1	150 x 190 cm	instrument table cover



Laparoscopy Bariatric Set

REF 1230611128**Number of packages:
1/5/10**


quantity	size	description
1	5 x 50 cm	op-tape
1	79 x 145 cm	Mayo stand cover, reinforced
1	150 x 190 cm	instrument table cover
1	240/215 x 350 cm	gynecological-abdominal drape with integrated sleeves



Extremity set


REF 120 311 / 1320611307
**Number of packages:
1/6/12**



quantity	size	description
2	9 x 50 cm	op-tape
1	79 x 145 cm	Mayo stand cover, reinforced
4	30,5 x 34 cm	hand towel
1	36 x 65 cm	stockinette
1	230 x 300 cm	extremity drape, elastic aperture 7 cm, absorbent patches 75 x 80 cm, tube holders
1	150 x 190 cm	instrument table cover



Hand/foot set


REF 120 325 / 1320611302
**Number of packages:
1/6/12**

quantity	size	description
1	9 x 50 cm	op-tape
1	79 x 145 cm	Mayo stand cover, reinforced
4	30,5 x 34 cm	hand towel
1	230 x 300 cm	extremity drape, elastic aperture 3,5 cm, absorbent patches 75 x 80 cm, tube holders
1	150 x 190 cm	instrument table cover





Hip set



REF 120 005 / 1320611338

Number of packages:
1/4/8

quantity	size	description
2	9 x 50 cm	op-tape
1	79 x 145 cm	Mayo stand cover, reinforced
1	36 x 120 cm	stockinette
1	100 x 150 cm	utility drape
1	150 x 300 cm	adhesive drape, absorbent patches 15 x 50 cm, tube holders
1	230 x 260 cm	adhesive split sheet, split 20 x 85 cm, absorbent patches 50 x 100 cm, tube holders
1	150 x 190 cm	instrument table cover




Knee arthroscopy set



REF 120 355 / 1320611345

Number of packages:
1/5/10

quantity	size	description
1	9 x 50 cm	op-tape
1	79 x 145 cm	Mayo stand cover, reinforced
2	30,5 x 34 cm	hand towel
1	36 x 65 cm	stockinette
1	230 x 300 cm	extremity drape, elastic aperture 7 cm, absorbent patches 50 x 100 cm, tube holders
1	15 x 240 cm	tube cover, perforated tip
1	150 x 190 cm	instrument table cover




Knee arthroscopy set



REF 120 035 / 132061B304

Number of packages:
1/4/8

quantity	size	description
2	9 x 50 cm	op-tape
1	79 x 145 cm	Mayo stand cover, reinforced
4	30,5 x 34 cm	hand towel
1	22 x 75 cm	elastic stockinette
1	150 x 190 cm	instrument table cover, folded
1	230 x 320 cm	knee arthroscopy drape, fluid collection pouch 105 x 90 cm with sieve and suction port, reinforcement wire and 2 elastic apertures Ø 5 cm
1	150 x 190 cm	instrument table cover




Abdomino-perineal set



REF 120 485 / 1320611512

Number of packages:
1/5/10

quantity	size	description
2	9 x 50 cm	op-tape
1	79 x 145 cm	Mayo stand cover, reinforced
2	30,5 x 34 cm	hand towel
1	215/240 x 350 cm	Gyn/Abdominal drape all-in-one with integrated leg covers and fluid collection pouch, with sieve and suction port, 2 adhesive apertures: top 25 x 30 cm, bottom 9 x 15 cm (with integrated cover drape if apertures not used)
1	150 x 190 cm	instrument table cover



Delivery set



REF 120 550 / 1320611509

Number of packages:
1/12

quantity	size	description
6	21 x 25 cm	hand towel
1	60 x 60 cm	underpad
1	60 x 90 cm	underpad
1	56 x 75 cm	baby sheet
1	75 x 120 cm	drape (stitching sheet)
1	90 x 150 cm	drape (for under the mother)




75 x 120 cm



60 x 90 cm




21 x 25 cm



56 x 75 cm




60 x 60 cm



90 x 150 cm



C-Section set I**REF 120 060 / 1320611506****Number of packages:
1/5/10**

quantity	size	description
1	9 x 50 cm	op-tape
1	79 x 145 cm	Mayo stand cover, reinforced
4	30,5 x 34 cm	hand towel
1	75 x 120 cm	baby sheet
1	180/250 x 330 cm	C-section drape, fluid collection pouch 72 x 72 cm, incision film, with aperture Ø 16 x 18 cm, 2 suction ports
1	150 x 190 cm	instrument table cover




Gyn/Cyst set**REF 120 055 / 1320611531****Number of packages:
1/5/20**

quantity	size	description
2	75 x 120 cm	leggings 75 x 120 cm
1	100 x 175 cm	gyn/cyst drape, adhesive aperture 9 x 15 cm
1	100 x 150 cm	instrument table cover

100 x 175 cm
Ø 9 x 15 cm


100 x 150 cm



75 x 120 cm

Vaginal set 1

REF 1230611535			Number of packages: 1/5/10
quantity	size	description	
1	75 x 90 cm	drape - 2-layers	
1	240/215 x 277 cm	drape gynecological with integrated sleeves	
1	150 x 190 cm	instrument table cover	




TUR-set



REF 120 072 / 1320611513


Number of packages:
1/6/12

quantity	size	description
1	9 x 50 cm	op-tape
1	2,5 x 30 cm	velcro tape
1	30,5 x 34 cm	hand towel
1	175/282 x 172 cm	TUR-drape all-in-one, finger cot, apertures 5 cm and 8 cm (adhesive), integrated fluid collection pouch 63 x 77 cm with sieve and suction port
1	150 x 190 cm	instrument table cover




Instrument Table Cover

REF	size	quantity
170 015 / 1320612401	100 x 150 cm	1/18/72
170 025 / 1320612402	100 x 190 cm	1/15/60




Mayo stand cover reinforced

REF	size	quantity
170 012 / 1320612404	79 x 145 cm	1/20/80




U-drape, absorbent patch, split 7 x 70 cm, 2-layers

REF	size	quantity
130 010 / 1320612101	150 x 250 cm	1/8/32





Drape - 2-layers

REF	size	quantity
140 010 / 1320612201	50 x 60 cm	1/45/180
140 015 / 1320612202	75 x 90 cm	1/25/100




Adhesive drapes - 2-layers

REF	size	quantity
1230612204	75 x 90 cm	1/100/200
1230612206	150 x 240 cm	1/20/40
1230612212	90 x 120 cm	1/65/130
1230612213	100 x 150 cm	1/40/80
1230612214	150 x 180 cm	1/55/55




Adhesive aperture drapes - 2-layers

REF	size	quantity
150 010 / 1320612304	75 x 80 cm, aperture 6 x 15 cm	1/25/100
150 015 / 1320612305	75 x 90 cm, aperture 6 x 8 cm	1/25/100
150 025 / 1320612306	150 x 180 cm, aperture 5 x 7 cm	1/14/56
150 030 / 1320612307	150 x 180 cm, aperture 6 x 15 cm	1/12/48
150 020 / 1320612314	120 x 150 cm, aperture 5 x 7 cm	1/20/80




Adhesive aperture drapes - 2-layers

REF	size	quantity
150 035 / 1320612308	230 x 250 cm, aperture 12 cm	1/7/28
150 039 / 1320612311	200 x 280 cm, aperture 15 cm	1/7/28
1230612320	120 x 150 cm, aperture 12 cm	1/40/80




Adhesive aperture drapes - 2-layers, variable

REF	size	quantity
150 013 / 1320612316	75 x 100 cm	1/20/80




Aperture drapes - 2-layers

REF	size	quantity
150 055 / 1320612301	50 x 60 cm, aperture 6 x 8 cm	1/45/180
150 065 / 1320612302	75 x 90 cm, aperture 6 x 8 cm	1/25/100




Ophthalmic drapes with incision foil

REF	size	quantity
1230618103	120 x 120 cm	1/35/140
1230618104	150 x 150 cm	1/20/80




120 x 120 cm



150 x 150 cm


Hand/foot drape

REF	size	quantity
130 020 / 1320612104	150 x 180 cm, aperture 3,5 cm	1/17/34




Extremity drape

REF	size	quantity
1230612103	230 x 300 cm, aperture 7 cm	1/10/20
1230612105	230 x 300 cm, aperture 3,5 cm	1/10/20




Tube covers perforated tip

REF	size	quantity
180 020 / 1320613202	15 x 240 cm	1/30/120




Dome Bags for X-ray

REF	size	quantity
180 055 / 1320613204	Ø 140 cm	1/10/40
180 050 / 1320613205	Ø 90 cm	1/18/72




Fluid collection pouches


REF	size	quantity
180 070 / 1320613106	30 x 40 cm, 1 chamber	1/55/220
180 072 / 1320613107	30 x 40 cm (30 cm x 15 cm + 25 cm), 2 chambers, suction and diathermy bag	1/55/220
180 073 / 1320613108	62 x 50 cm, sieve and outlet	1/11/44



30 x 40 cm



30 x 40 cm



62 x 50 cm

Hand towel

REF	size	quantity
1230613111	29 x 30 cm	1/90/720



Op-tape

REF	size	quantity
1230613112	9 x 50 m	1/120/480




Armrest cover

REF	size	quantity
1230613124	20 x 55 cm	1/80/320




Arm cover

REF	size	quantity
1230613110		1/100/200




Incision film

REF	size	quantity
1230610101	15 x 20 cm	1/120/480
1230610102	30 x 20 cm	1/100/200
1230610103	45 x 20 cm	1/15/240
1230610104	40 x 35 cm	1/40/40
1230610105	45 x 50 cm	1/30/30
1230610106	45 x 65 cm	1/30/30




Urinary Catheterization Set 672 - Sterile

REF 1230111672			Package quantity: 1/84
Quantity	Size /	Description	
4	Ø 30 mm	Non-woven ball	
2	5 x 5 cm	Non-woven swab, 4-ply	
1	13 cm	Plastic forceps, blue	
1	10 ml	Syringe BD PosiFlush SP, 0,9% NaCl	
2	Size M	Nitrile examination gloves, non-powdered	
1	50 x 45 cm	Utility drape	




Suture Set 356A sterile**REF 1230111356A****Package quantity:
1/23/46**

Quantity	Size	Description
1 pc	45 x 75 cm	Hospidrape
1 pc	13,5 x 27 x 25 cm	Bowl, blue
1 pc	50 x 60 cm	Adhesive drape with aperture 6 x 8 cm
5 pcs	7,5 x 7,5 cm (4-ply)	Nonwoven swab
1 pc	14 cm	Needle holder
1 pc	11 cm	Scissors Iris, sharp
1 pc	12 cm	Forceps Adson



Dressing Set 160V sterile**REF 1230111160V****Package quantity:
1/80/160**

Quantity	Size	Description
5 pcs	20 / 19 cm	Gauze ball
1 pc	13 cm	Plastic forceps, green
5 pcs	7,5 x 7,5 cm (8-ply)	Gauze swab



Central Venous Catheterization Set 164V sterile

REF 1230111164V
**Package quantity:
1/20/20**

Quantity	Size	Description
2 pcs	Size M	Latex gloves, powder-free
1 pc	75 x 90 cm	Hospidrape
1 pc	10 x 12 cm	PU transparent tape Dermafoil
1 pc	21 x 14 x 4,5 cm	Bowl, transparent
2 pcs	5 x 5 cm (8-ply)	Gauze swab
6 pcs	7,5 x 7,5 cm (8-ply)	Gauze swab
5 pcs	30 / 30 cm	Gauze ball
1 pc	10 ml	Syringe, 2-piece, Luer
1 pc	20 ml	Syringe, 2-piece, Luer
1 pc	0,7 x 30 mm, 22G	Needle, black
1 pc	1,2 x 40 mm, 18G	Needle, pink
1 pc	13 cm	Kelly clamp, straight
1 pc	13 cm	Plastic forceps, green
1 pc	14 cm	Pean clamp, straight
1 pc	11 cm	Scissors Iris, straight
1 pc	No. 11	Disposable scalpel
1 pc	75 x 90 cm	Adhesive drape with aperture Ø 6 cm, with perforation

